Ritter ProCure

C20-0058

2201

Polymerisationslicht

Curing light

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- Instruction for use -

Ritter® Concept GmbH
Freiburger Str. 45
88400 Biberach, Germany
Fon: +49 7351 52 925-10
Fax: +49 7351 52 925-11

e-mail: info@ritterconcept.com
Internet: <http://www.ritterconcept.com>

 **Ritter**®
[THE DENTAL EXPERTS]

Packliste

Inhalt

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Packing list

Content

- Handpiece with LED 1 set
- Optical fiber 1 pcs
- Light hood 1 pcs
- Instruction for use 1 pcs

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| Icon | Legend |
|---|---------------------------------|
|  | Follow the instructions for use |
|  | Type B applied part |
|  | Protection class II |

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Intended use

This product is used to restore teeth and solidify material for whitening teeth.

Principle

Ritter ProCure (build-in) adopts the principle of ray radiation to solidify light-sensitive resins in a short amount of time.

Structure and components

Ritter ProCure (built-in) is composed mainly of high power LED, optical fiber, main unit and power supply line.

Technical Specifications

| | |
|------------------------|--|
| Type: | Curing light |
| Rated voltage: | 24 V ~ |
| Rated input: | ≤8 W |
| Electrical protection: | Class II |
| Protection class: | Type „B“ ⚡ applied part |
| IP-Code: | IPX0 (no protection against foreign objects and water) |
| Light source: | blue light (385-515nm) |
| Power: | from 1000 mW/cm ² to 1200 mW/cm ² |

Working conditions:

| | |
|--------------------------|-------------|
| Environment temperature: | 5°C to 40°C |
| Relative humidity: | 30% to 75% |

| | |
|------------|--------------|
| Dimension: | 26x25x260 mm |
| Weight: | 135 g |

Attention

Do not use in area of flammable anesthetic mixture with air, oxygen or nitrous oxide.

Technical modifications reserved!

Date of issue: 2022-03-30
Original language: German

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Install and uninstal

1. Connect the LED power supply line with the power (15V-24V AC) of your dental unit. Thighten the nylon thread to the fixation point of your dental unit, then it is ready for operation.

Note:

When installing the LED, be sure the power is cut off. The two power wires should be a little longer than the nylon thread to keep the power wire safe.

2. Take the red protection cap off the fiber optic and then insert metal cap of the fiber optic into the curing light's installation hole. (You must insert the parts on the end point by twirling, be sure not to insert on an angle)
3. Install the glare protection on the top of the fiber optic.
4. To uninstal the LED, just reverse the procedure above.

Operation

1. You can choose one operating mode by pressing the mode switch button (Pay attention to the LED display on the handpiece)
 - (1) **Full power:**
The curing light radiates at max. intensity
→ light duration: 10s
 - (2) **Ramping:**
The curing light intensity increases continually (after 5s it reaches the max. intensity).
→ light duration total: 10s
 - (3) **Pulse:**
The curing light works in pulse mode.
→ light duration total: 10s

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Note:

The light intensity of this (LED driven) curing light is far higher than a halogen lamp driven curing light.

The depth of the composites cured by the light in 10 seconds is at least 4 mm.

4. The optical fiber of the curing light can be removed and autoclaved (at 134°C / 0.22MPa).

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Note:

- To avoid a wrong radiation position or improper curing at clinical use, the curing light needs to be set vertical on composite materials/resins.

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| | |
|--|--|
| | Do not light the (blue) curing light into the eyes! |
|--|--|

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Daily maintenance

- After every operation, shut off the electricity and check the fiber optic whether there is any residual resin, wipe off if any.
- Wipe the surface with drinking water or distilled water.

Note:

- **Do not immerse in any liquids.**
- **Only the fiber optic can be autoclaved (in steam sterilizer).**

Reprocessing instruction for Ritter ProCure to be performed after each patient**Warning**

- Sterilization must be performed even when using the optical fiber for the first time
- Inappropriate sterilization is dangerous for patients and operators
- Do not perform spray sterilization
- Do not perform cold chemical or hot air sterilization
- The appointed staff must be specialized and trained
- Use the disinfectant in accordance with the specifications set by the manufacturer indicated on the label/ instructions for use
- Do not use chlorine-based solutions
- When sterilizing in an autoclave, do not exceed the maximum permissible filling quantity

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Troubleshooting

| Description | Reason | Troubleshooting |
|------------------------------|--|--|
| No indication, no activity | The LED is not connected with the power or the power switch is off. | <ol style="list-style-type: none"> 1. Check the connection of the LED and the power. 2. Set power switch „ON“. |
| Light intensity insufficient | <ol style="list-style-type: none"> 1. Optical fiber isn't inserted in base position. 2. The optical fiber itself has a crack. 3. The resin remains on the surface of the optical fiber. | <ol style="list-style-type: none"> 1. Please install the optical fiber again. 2. Please change the optical fiber. 3. Please clean the optical fiber - remove residual resin |

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Warranty

1. We offer one year repair to those products that can't work because of quality problems, counting from the sales date.
2. There is no warranty claim for damage caused by improper use.

For any questions please contact us via e-mail: sales@ritterconcept.com.

Storage and transportation

- This device should be handled carefully.
- Do not shake, protect against shocks (do not drop).
- Installed or stored at shadow, dry, cool and ventilated places.
- Don't store together with substances that are flammable, poisonous, caustic or explosive.

Transport conditions:

Environment temperature: -20°C to +55°C

Relative humidity: 10% to 93%

Atmospheric pressure: 700 hPa to 1060 hPa

- Keep it away from sun or rain/snow during transportation.
- Avoid shocks.
- Don't mix it with dangerous substances during transportation.

Disposal

There are no harmful elements in this device. It can be disposed according to the local laws.

Note

The illustrations in this manual may differ from the actual product.