Ultimate Comfort

C20-0087
2302

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Dental work place
(electric and hydraulic version)

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Dear user!

Our congratulations on the choice of your new Ritter product and many thanks for your confidence! This modern device was produced by our development- and production engineers with great carefulness.

You and your team will gain absolute joy and benefit of your new Ritter device, when following the tips and advices, given to you during the handing over demonstration.

Despite of a careful introduction into the operation, smaller or greater questions will arise in the course of time. This operating manual shall give you answers to these questions. To make the working easily for you and to be sure, that all technical or legal regulations and the maintenance notes receive their attention, we recommend a careful examination of this manual for you and your responsible employees too.

If you face technical problems, our service hotline
+49 (0) 7351 52925-37 is available for you.
And now a lot of joy with your Ritter equipment.

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Date of issue: 2023-10-11
Original instruction for use

Version: C20-0087-2302

Explanation:

C20 = Code for user manual
0087 = Document identification
2301 = Version number consisting of:
(20)23 current year and
02 current revision of document


Abb.: Dentale Behandlungseinheit "Ultimate Comfort"
Fig.: Dental treatment unit "Ultimate Comfort"

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No.	Description
1	Assistant element
2	Spittoon bowl
3	Cup filler
4	Operating lamp
5	Table arm
6	Dentist element
7	Dentist keyboard
8	Multifunctional foot controller
9	Patient chair

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Patient chair

Patient chair - dimensions.....	42
Control elements, manual control.....	43
Program control.....	44
Head rest.....	46
Operating lamp – SunLite.....	47
Operating lamp – CoolLite.....	49

Maintenance

Daily.....	52
Troubleshooting.....	58

Appendix

Setting up drawing.....	62
Air-water-schematic.....	63
Wiring diagram – UC electric version.....	64
Safety inspection.....	66
EMC-declaration.....	69
RiTron.....	72
RiTren.....	74
Terms of Guarantee.....	75

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Dental work place

Type:	Ultimate Comfort
Version:	electric
Rated voltage:	230 V 50/60 Hz
Rated input:	900 VA
Fuse	
in device installation:	T 10 AL 250 V
Protection class:	I
Type applied part:	B
IP-Code:	IPx0
Air source:	
Air pressure:	>5.5 bar (>550 kPa)
Air flow:	50 l/min
Air filter:	20µm
Air source requirement:	oilfree, dry, without bacteria
Water supply:	
Water pressure:	2 ... 4 bar (200...400 kPa)
Flow:	> 10 l/min
max. waterconsumption:	2.8 to 4.5 l/min

Patient chair (electric version)

Mode of operation:	DAB* (2min/18min)
(Intermittent Duty Time Rating):	On: 2min Off: 18min)
Kind of movement:	electric system
Maximum height (seat cushion ⇄ floor):	800 mm
Minimum height (seat cushion ⇄ floor):	450 mm
Maximum stroke:	350 mm
Backward angle of backrest:	115°-180°
Up and down sliding stroke of head rest:	200 mm
Maximum load:	150 kg

Water unit

rotational angle of spittoon bowl:	≥ 180°
Water filter:	75 µm
Heating	
Rated voltage:	24 V 50/60 Hz
Rated input:	120 VA

Dental operating light „SunLite“

Typ:	LED
Rated voltage:	17-24 V ~ 50/60 Hz
Rated input:	20 VA
Mode of operation:	DAB**
Protection class:	II
Illuminance	8000-35000 lx
Colour temperature:	4200-6000 K
Colour rendering index (CRI):	>85
Lifetime cycle:	50000 h

Dentist element


Maximum loads:	3 kg
----------------	------

DAB* = continuous operation with intermittent load.
The permissible load times are adapted to the working mode of the dentist (please refer also to the operating manual for the individual units and/or instruments).

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Explanations to name plates

Icon	Explanations to type plates
	Serial number
	Medical Device
	<u>CE-Approval</u> according to Directive 93/42/EEC for medical products
	The symbol of the crossed waste container on wheels indicates, that Ritter products are electrical equipment according to Directive 2012/19/EU. Within the EU-Countries these products have to be separately disposed. Ritter wants to encourage you, to minimize possible effects of waste to environment and health outside of the EU, too. Please follow the local directives for disposal of waste and make use of the separate collection of electrical equipment.
	There are specific warnings or precautions related to this product.
	Type B applied part
	Follow the instructions for use
	Manufacturer/Date of manufacturing (according ISO 15223)

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Operating conditions


Temperature: +10°C to +40°C
rel. humidity: 30% to 75%
Barometric pressure: 700 hPa to 1060 hPa

Transport and storage conditions

(original packing)
Temperature: -20°C to +55°C
rel. humidity: max. 93%
Barometric pressure: 500 hPa to 1060 hPa

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Intended purpose

This dental treatment unit is used for diagnosis, therapy and dental treatment in humans. It positions the patient, illuminates the oral cavity and provides the drive for the dental instruments.

Indication

Not further specified.

Contraindication

Patients who cannot be positioned on the chair and cannot maintain the position required during treatment.

Patient target group

There are no specific restrictions based on age, weight, health status or ethnicity as long as patients are able to follow the user's instructions. The decision as to whether a patient is eligible for treatment with the product rests solely with the intended user. The maximum load of the patient chair must be observed.

Intended user

Trained and instructed dental professionals.

Complication

- Risk of crushing: When moving the backrest, the patient's arms must not be behind the back or moving dentist or assistant element
- Reversible glare if the LED treatment light is handled incorrectly.
- Risk of infection due to transmission of pathogens through cross infection
- Risk of infection due to transmission of pathogens
- Risk of overloading the dentist element, assistant element or patient chair

Clinical benefit

Both patient and user benefit from the ergonomically arranged components of the patient chair, dentist element, work field light, water unit and assistant unit.

The patient is positioned by the patient chair so that access to the oral cavity is possible for the user. The patient chair is adjusted in height and the inclination of the backrest is set accordingly.

Dental treatment is performed with additional instruments and devices. These are supplied with electrical voltage, compressed air, vacuum and water by the dental treatment unit.

Serious incidents

All serious incidents related to the product must be reported to the manufacturer and the competent authority of the Member State where the user is established.

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Our equipment complies with the requirements of of Medical Devices Implementation Act (MPDG) and Directive 93/42/EEC for Medical Devices, revised by EU-Directive 2007/47/EG from 2010-03-21.

It complies with the VDE-requirements for electromedical devices DIN EN 60601, part 1 / DIN VDE 0750, part 1, too. It is suited for connection to supply lines installed according to the DIN VDE 0100-710 regulations for medically used rooms.

	This unit is not explosion-proof.
	This dental treatment unit must not be modified or opened without permission of the manufacturer.
	We can only accept responsibility for the safety properties of our equipment, if installation, maintenance, repair and modifications are carried out by our workshop or by a person explicitly authorized by us. Further, the equipment has to be correctly handled according to our operating manual.
	In case of failure, structural parts, that can affect the security of the device, are to be replaced by original parts only.
	In case of repair, the technician must issue a certificate as to the kind and extent of work performed and, if applicable, indicate any changes of nominal data and operating range. The certificate must contain the date of repair, the name of the executing company and a signature (according to MPDG, DIN VDE 0750 / DIN VDE 0751).
	Close to the equipment do not operate any devices producing electro-magnetic emissions.
	An ergonomic way of working must be chosen in order to prevent musculoskeletal disorders.
	Before using for the first time, check the attachment of the ground anchoring!

Safety-relevant inspections

The safety standards for medical products and all legal obligations are to be fully observed by the user.

Scope and schedules of safety-relevant inspections are as follows:

Every two years the following safety-relevant inspections according to DIN EN ISO 62353 have to be performed: see Form „Safty inspections – Dental treatment unit“.

We recommend to keep an equipment control log for documenting the results of the safety-relevant inspections.

The unit lifetime cycle is: 10 years.

	Multimedia For components like Intraoral camera, monitor, note the separate user manuals for these devices!
---	--

Disposal

According to law, electric equipment is to be disposed by authorized disposal firms only.

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Please always follow the following safety precautions to avoid injuries to persons or damages to the unit:

	WARNING!
	Always shut off the water supply, air supply and the main power before performing any cleaning or maintenance.
	Install the dental unit on level ground.
	Do not install the unit near radiators or heaters.
	Do not place or hang any heavy things on the unit.
	Do not connect unit with incompatible wires which may affect its performance, or even cause fire or electric shock.
	Do not place anything in the unit so as to avoid electric shock or fire. Do not spill any liquid into the inner part of the unit.
	Do not disassemble the unit so as to avoid electric shock. Contact technical service personnel whenever the unit needs repairing.
	WARNING: To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth
	Do not open or remove the cover plate so as to avoid injuries caused by dangerous voltage or other dangers. An improper assembly may cause electric shock during the subsequent operation.
	Please switch off the power and seek help from technical service personnel if: (1) any break-down or wear-down of the plug or cable occurs. (2) unknown liquid spills into the inner unit or water flows into the inner circuits. (3) the unit malfunctions, even if it is operated according to the instructions.
	If possible, turn off the power supply during thunderstorms.
	Please never exceed this units weight limit (maximum load of the dental seat: 150kg).
	The maximum load of the table tray attached to the dentist table is 3.0 kg, which should never be exceeded.
	The manufacturer shall not be responsible for any damage or injury caused by mishandling, neglect or improper use.
Before leaving the surgery, shut off the water supply, air supply and the main power switch.	


Abb.: Arbeitsbereich "Ultimate Comfort"
Fig.: Working area of "Ultimate Comfort"


Abb.: Arbeitsbereich "Ultimate Comfort - elektrische Version"
Fig.: Working area of "Ultimate Comfort - electric version"


Abb.: Arbeitsbereich "Ultimate Comfort - hydraulische Version"
Fig.: Working area of "Ultimate Comfort - hydraulic version"

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Translation of user manual

When imported into countries whose native language is not English or German, the operating manual are translated into local languages. In this translation, replace the words „original operating manual“ by the following notice: „translation of the original operating manual.“

The Ritter® Concept GmbH is not directly or indirectly liable for mistakes or failures in this documentation as well for the given benefits of the use with this operating instruction!

The illustrations in this manual may be able to differ from the actual unit in your system.

The duplicating of this documentation without written permission of Ritter® Concept GmbH is prohibited!

Accompanying documents are written at a level consistent with education, training and other needs of individuals for whom they are intended.

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- Grundplatte (Holz)
- Kartonverpackung
- PE-Kunststoffhaube

Recycling

All packing material does not harm the environment and can be completely recycled.

- wooden base plate
- carton box
- polyethylene bag

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We recommend you to read this operating manual very carefully in order to have full knowledge of all the advantages this dental workplace is offering you.

The essential attributes of all models are:

- 5 outlets for instruments (air/water)
- cut off air/water switch for safety during changing drills

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
Model "H" (hanging hoses)

- instrument carrier with hanging hoses smooth running swivel arm
- pneumatic arm brake
- integrated tray
- photoelectric sensor in instrument holder

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Model "S" (swinging hoses)

- Instrument carrier with swinging lever
- smooth running swivelarm
- pneumatic brake

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
Abb.: Netzschalter
Fig.: power switch

Startup

The unit is activated via the power switch (1) at the patient chair pedestal. After switching-on the control indicator of the power switch lights.

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It is recommendable to actuate the hand pieces, the cup filler and the spittoon bowl rinsing before the first daily treatment is done in order to rinse out stale water.

	<p>For safety reasons, switch off the main switch when leaving the surgery!</p> <p>Therefore press main switch (1). After switching-off, the control indicator of the power switch turns off.</p> <p>Using stop-cocks to shut off air and water inlet.</p>
--	---

If suddenly power outage switch off the main switch (1) and take all handpieces in their instrument holder. When electrical power is turned on again, activate the unit via the power switch (1).

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Arm brake

Moving the dentist element:

- touch the dentist element on the handle (2)
- press the button (1) and move dentist element in the desired position
- release button (1)

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Control elements patient chair

A	Pedal / Foot lever
If all instruments (except the syringe) are in their holders.	
2	Slide tip switch to the <u>left</u> → Backrest DOWN
	Slide tip switch to the <u>right</u> → Backrest UP
	Slide tip switch <u>backwards</u> → Patient chair UP
	Slide tip switch <u>forwards</u> → Patient chair DOWN
3	press key <u>shortly</u> : Automatic movement to „Spitting position“, actua- ting again moves chair into last treatment position
	press key longer than <u>2 seconds</u> an acoustic signal sounds, Automatic movement to programmed „Position A“
4	press key <u>shortly</u> : Automatic movement to „Exit position“
	press key longer than <u>2 seconds</u> an acoustic signal sounds, automatic movement to programmed „Position B“

Control elements for instruments

3	press key <u>shortly</u> : activates the spray water at instruments
4	with <u>pressed</u> key activates the chip blower-function
5	push lever to the <u>right</u> : starting the activated instrument and increase/decrease the speed/power
6	maximum speed/power of the instrument
7	Press the pedal. Start the active instrument and increase/decrease the speed/power

Programming the patient chair positions:
→ see page 36.

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
Functions:

Water, air or spray are switched by hand. Supply of the media lasts as long as the corresponding key of the hand piece is pressed. Either air or water can be obtained by actuating the corresponding key. Actuating both keys at the same time will supply spray.

Operation of 3F syringe:

- 1 Key for water delivery
- 2 Key for air delivery
- 1+2 Delivery of spray

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1	Key for water delivery
2	Key for air delivery
3	Spigot ring
4	Cannula
5	Three-way Syringe
6	hand piece

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 Removable and autoclavable cannula (4).
(sterilizable in steam sterilizer at 134°C)

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Spittoon bowl rinsing

The spittoon bowl rinsing starts when pressing key „Spittoon bowl rinsing“ and stops automatically after a programmed time.

Premature stop during bowl rinsing is possible by pressing the key once more.

Programming of spittoon bowl rinsing time:

- press key  longer than 5 seconds
→ the indicator light blinks
- press key  shortly
- when the desired rinse time is reached, press key  again to finish the spittoon bowl rinsing
→ indicator light turn off and the programming mode is finished

Cup filler

The cup filling starts when pressing the key “Cup filler” and stops automatically after a programmed time.

Premature stop of cup filling can be released by pressing key “Cup filler” once more.

The recommended filling level is two thirds of the glass or plastic cup. Overflowing water will run off through the cup support plate to the drain.

Programming of cup filling time

- press key  longer than 5 seconds
→ the indicator light blinks
- press key  shortly
- when the desired filling level is reached, press key  again
→ indicator light turn off and the programming mode is finished

Cup fill heating

- press key  shortly
→ the indicator light blinks
When the water reaches the preset temperature the indicator light lights continuously.

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X-Ray viewer

- press key  shortly
- after 3 minutes, the x-ray viewer turns off automatically
- Luminance $\approx 1600 \text{ cd/m}^2$

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Instrument light

- take out the instrument from the holder
- press key  shortly to switch the light at the handpiece ON/OFF

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Operating lamp

- press key „Operating lamp“ to switch the operating lamp ON/OFF

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Scaler

- screw the required scaler tip into the hand piece (clockwise)
- insert the scaler hand piece into the hose connector
- start the scaler by activating the foot control

Adjust the power of the scaler

- rotate the power regulator (2) at bottom of the dentist element:
 - decrease power: rotate clockwise
 - increase power: rotate anti-clockwise

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
Adjust water flow for scaler

- by using the corresponding water regulator (1) underneath of the dentist element
 - see page 26

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


For each instrument (except the syringe) there is the possibility to regulate the spray water quantity by hand.

The regulators can be found at the bottom of the dentist element (1).

The order of regulators depends on the order of the instruments on the tableau.

In order to regulate the water, take the according instrument from its holder and actuate it with spray water function. Now you can adjust the water flow by turning the respective regulator as illustrated in the picture below.

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
Air and Water pressure adjustment of handpieces:

There are air and water regulating valves used for adjusting the water vapor.

1	Luftdruckregler für 3. Bohrinstrument	air pressure regulator for 3. th burr instrument
2	Luftdruckregler für 2. Bohrinstrument	air pressure regulator for 2. nd burr instrument
3	Luftdruckregler für 1. Bohrinstrument	air pressure regulator for 1. st burr instrument

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
Adjustment drive air for handpiece:

Low speed pneumatic handpiece

No-load speed

- shall be no lower than 14000 min⁻¹ if the air pressure is 2.5 bar (250 kPa, driving air)

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High speed turbine

No-load speed

- shall be no lower than 30000 min⁻¹ if the air pressure is 250 kPa (driving air)
- a air pressure regulator is at the bottom of the dentist element
→ reduce the air pressure by rotating the regulator clockwise (1, 2, 3)

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Enable/disable the prompt tone function

- press key  longer than 5 seconds to enable/disable the prompt tone function
→ a buzzer beep sounds
→ the prompt tone function of the key is enabled/disabled

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1	Big suction
2	Curing light
3	3-function syringe
4	Small suction

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Big suction (1)

- take out the suction handpiece
→ suction starts automatically
- put the suction handpiece back into the instrument holder
→ suction will stop automatically after 3 seconds

Small suction (4) (small suction)

- take out the suction handpiece
→ suction starts automatically
- put the suction handpiece back into the instrument holder
→ suction will stop automatically after 3 seconds

 **Suction hose and suction handpieces are expendable materials.**

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 **Bottle can be unscrewed only when it is depressurized.**

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1	air break switch - bottle compressed air <ul style="list-style-type: none"> • „Bottle OFF“ → Bottle in pressure-free mode • „Bottle ON“ → Bottle is now under pressure
2	Toggle - water source <ul style="list-style-type: none"> • „Tap water“ → water from water pipeline • „Bottle“ → water from bottle

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The use of distilled water is recommended.
If there are impurities in tap water, it is necessary to use the distilled water from bottle.

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General directives

Manufactured under observance of: EN 13485, Medical Devices Implementation Act and the directive 93/42 EEC.

This supply unit corresponds to the norm DIN 1988-100 / DIN EN 1717 for protection of drinking water to prevent a back flow of the process water (of peripheral devices and dental workplaces) to the public water net.

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	<p>IMPORTANT: For safety reasons, the service life of the bottle included in this system is limited to a period of six months. After this time, it is absolutely necessary to replace the bottle with a new one. Otherwise, the bottle may be destroyed under pressure, resulting in damage to health and/or damage to your treatment center.</p>
--	--

Clues for disinfection of the dental unit

With help of intended use of the Ritter® AquaBottleSystem the aquiferous systems of the Dental unit can be disinfected periodically. The water disinfection must be performed with the detergent/disinfectant names RiTron. Otherwise a proper disinfection can not be guaranteed.

The following points should be noted:

- Because of the device compatibility use only the specified disinfectant RiTron.

	<p>IMPORTANT: After every successful disinfection performing RiTron must be rinsed out of the system <u>completely</u>. Scour the water system so that the RiTron will be rinsed out of the water unit. That is indicated by a loss of the blue coloration of the rinsing water. When clear water comes out from the instrument hoses and from the Cup filler, the wash cycle may be terminated – not earlier!</p>
--	---

- For an optimal disinfection the disinfectant must stay 24 h or longer in the water system (e.g. during weekend).
- The blue RiTron may only be poured into the containers provided.
- It must be ensured, that there is enough RiTron in the bottle before the cleaning/disinfection start. In the case of premature emptying air would be sucked in and the system would not be filled with disinfectant completely. The system would not be filled with disinfectant completely. Only a water unit, which is filled with RiTron completely can be disinfection dependable! Furthermore RiTron can be splashed out uncontrolled squirting because of air in the water hoses (for example when activating the rinse glass filler)
- Disinfection must be performed according to hygiene regulations in dental practices (using protective gloves, protective glasses and so on).

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Filling of the bottle

Follow the instructions ‚change water bottle‘ on page 32. Fill the bottle with the correct amount of RiTron from the storage container and handtightly rescrew the bottle to the receptacle.

When the bottle is switched on, a pressure builds up. This is clearly audible due to the incoming air.

→ The dental treatment unit is now ready for a disinfecting process using RiTron.

Performing the disinfection

RiTron should be filled into the water-bearing systems on Friday or before lasting several days of work break. Ensure that enough RiTron is in the bottle before beginning any disinfecting process. For one disinfecting process 300-600ml RiTron is needed (**Observe the fill level marking on the bottle!**). Every water-bearing system must be rinsed as long as the blue disinfecting agent RiTron is leaking out.

- put all instruments on maximum flow
- Detach the turbine, the handpieces and the contra-angle handpieces (adhere to the instrument manufacturer's user manual)
- Start with the cup filling for about two - three times
- Activate the single instruments via hand or foot controller with active water
- If a rinsing system is connected to the suction hoses or to the separation automatic, activate every suction hose for about 60 s
- External connected devices (e. g. Airflow, if existent at the unit) are also to be rinsed!
- Every external device must basically be rinsed! (A therefore necessary plug with hose for the rinsing of the external device connection can be ordered from the manufacturer or the dental depot)
- When the disinfection procedure has been finished, switch off the dental unit and set the switch at the bottle support to ‚Bottle OFF‘ position. The bottle will be depressurized and can be filled with fresh drinking water and 1% of RiTron before starting to work

Dosing specification

Bottle volume: 1.0 l ⇔ add 10 ml RiTron from the measuring cup.

When working with RiTron, wear rubber/PVC gloves and tight fitting eye protection.

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	<p>IMPORTANT: <u>All</u> connected water hoses must be rinsed in every case! Even if they were not used in the week before. Remaining water is the ideal precondition for increasing of germs!</p>
--	---

Placing into operation of the dental unit

The blue disinfecting agent has to be rinsed out on Mondays or after several days lasting breaks of work. Every water-bearing system must be rinsed as long as clear water leaks out.

The following steps have to be performed:

- Fill the bottle with fresh drinking water and 1% RiTron. Set the switch at the bottle support to **On** position.
- Switch-on the dental unit.
- Activate the spittoon bowl rinsing and rinse the spittoon bowl for about 2 min to rinse out the stagnating water in the house main water line.
- Actuate the cup filler for about 8 times.
It is useful to empty the 8 cup fillings with the suction system to check its function.

	<p>Ensure that the handpieces and contra-angle handpieces are detached before starting to rinse. This prevents that the nozzles are being blocked with residues.</p>
--	---

- Rinse through every single instrument hose for about 2 min.
- Externally connected devices and external water connections have to be rinsed, too!
- When clear water is leaking out of the hoses, the turbines, hand- and angle pieces can be connected. When connecting external devices, make sure that the connectors have been disinfected to prevent the entry of new germ by the quick couplings.
- Fill the bottle again with fresh drinking water.
For Disinfection of the 1.0 l water in the bottle add 10 ml RiTron (≈ 1 dosing cup).
- If foaming appears on the turbine spray during treatment, the RiTron dosage must be reduced to up to 5 ml RiTron per 1.0 l water, depending on the water quality.

The dental unit is now ready for operation.

The periodic weekly disinfection should be registered into the included hygiene plan for self control.

 **The water outlet of the cup filler has to be cleaned in periodic intervals.**

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
Error Source

No pressure is establishing within the bottle and no disinfection agent leaks out of the instrument hoses!

- Check the correct placing of the bottle at the bottle support.
- The pressure lines have a leak → Inform a technician.
- Set the switch to „Bottle OFF“ position, tighten the bottle a little bit more to bottle support and switch-on again. (Attention: Don't tighten the bottle too firmly!)

RiTron was spilled out during the disinfection process!

- Clean the dirty area with available disinfection agents.
- It has to be noticed that the dirty objects/areas are impassible again the uses disinfection agent.

During the rinsing with disinfection agent, only little amounts of RiTron leaks out of the instruments!

- The air pressure of the house connection is too less so that inside the bottle not enough pressure is built up (Notice the limit value!).

During the rinsing with disinfection agent, only little amounts of clear water leaks out of the instruments!

- It was forgotten to detach the instruments from the hoses before starting to rinse.
 Residues (biofilm) of the hoses are blocking the instruments.

Spare parts:

Component/Spare part	Order No.
Gasket pressure bottle UC	F000034540
Replacement bottle 1.0 l	K001123026

Consumables

Description	Order No.
RiTren: 5 l Container (One 5 l container is included in the scope of delivery!)	U000003910
RiTren: 1 l bottle with dosing cup (one bottle is included in scope of delivery!)	U000003911

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Important safety regulations

Please carefully read the notes!

- | | |
|--|---|
|  | The bottle is pressurized! (ca. 2.0-3.8 bar, depends on the type of device) |
|  | Don't manipulate the bottle with sharp objects! |
|  | Use only original bottles from Ritter® Concept GmbH. |
|  | Before opening the bottle is necessary to check that the switch (valve lifter) is set to Off (see page 32). |
|  | After filling the bottle with water always screw the bottle in the correct way into the bottle holder! |
|  | The bottle must be screwed only hand-tight to the bottle support. |
|  | After disinfection clean the thread of the bottle and bottle holder regularly. |
|  | If the bottle has any mechanical defect, it has to be replaced by a new one. |
|  | For safety reasons, the use of the bottle is limited to a period of 6 months . After that the old pressure bottle has to be replaced by a new one. |

 In non-compliance with these safety regulations the company Firma Ritter® Concept GmbH shall not be liable for any resulting damages!

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Instructions for sample taking

- The samples must be extracted separately at each instrument. A mixing of the water samples does not give information about which instrument is the causing the exceeding of the limit values.
- Before sample taking all hoses (including cup filler) have to be rinsed for about 2 min to fill the hoses with fresh water.
- The water leak out on positions of the endings of the hoses (with or without plugged hand and angle pieces) and the cup filler has to be disinfected by wiping. Let the disinfectant soak long enough - as prescribed.
- Extract the necessary amount of water and transport it as fast as possible to according place (approved laboratory).

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Program control:

Treatment position A

- is called by key „Treatment position A“ at the dentist keyboard or foot controller
- the operating light will be turned on simultaneously

Treatment position B

- is called by key „Treatment position B“ at the dentist keyboard or foot controller
- the operating light will be turned on simultaneously

Spitting position/last position

- press key „Spitting position“ at the dentist keyboard, assistant keyboard or shortly press left tip switch of foot controller → automatic movement of patient chair to the „Spitting position“, spittoon bowl swivels inwards and the spittoon bowl rinsing starts; the operating lamp switches off
- actuating again moves chair into last treatment position, spittoon bowl swivels outwards; the operating lamp switches on again

Exit position

is called by key „Exit position“ at the dentist keyboard, assistant keyboard or by pressing shortly the right key at the foot controller

Emergency position

is called by key „Emergency position“ at the dentist or assistant keyboard (automatic movement of patient chair to the „Emergency position“ and the operating lamp switches on)

Programming the patient chair positions:

Treatment position A

- move the patient chair manually into the desired position
- actuate and keep actuated the program key  on the dentist or assistant keyboard longer than 5 s till a single beep sounds

Treatment position B

- move the patient chair manually into the desired position
- actuate and keep actuated the program key  on the dentist or assistant keyboard longer than 5 s till a single beep sounds

Spitting position

- move the patient chair manually into the desired position
- actuate and keep actuated the program key  on the dentist or assistant keyboard longer than 5 s
- a single beep sounds and the current position of patient chair will be stored

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Emergency position

- move the patient chair manually into the desired position
- actuate and keep actuated the program key  on the dentist or assistant keyboard longer than 5 s
- a single beep sounds and the current position of patient chair will be stored

Exit position

- move the patient chair manually into the desired position
- actuate and keep actuated the program key  on the dentist or assistant keyboard longer than 5 s till a single beep sounds

Reset to factory settings

- press key  longer than 5 seconds
 - the patient chair will move to the highest position, then to the lowest position
 - all 5 program memories for patient chair positions are reset automatically
 - a signal sounds - the reset to factory settings is finished

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Head rest

Adjust the headrest:

The headrest can be adjusted through the double joint using the rotary knob(1):

- turn the rotary knob (1) counterclockwise to loosen the joints
- adjust the headrest to the required position
- turn the rotary knob (1) clockwise to lock the joints of the headrest

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Right Handle

Move armrest UP

- during the treatment, move the end of the armrest handle upward
- a click sound indicates that the handle has reached the fixed treatment position

Move armrest DOWN

- move the end of the armrest handle upward to the highest position
- then the handle is movable to the lowest position

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Warning

During operation process, no tangible items shall be placed within the operational stroke of the unit!

- When the patient chair or backrest move down, the backrest will automatically stop as for obstruct encountering and moves up for a short distance.
- When the patient chair moves, press any key to stop the movement.
- When the dental handpiece stops working, the spray will delay about 2s and form automatic back-suction resistant function to effectively prevent cross infection.
- If a handpiece is removed from the instrument holder and in operation, the movement of the patient chair is blocked

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Operating light is fixed on flexarm mounted to the treatment unit.

The special design of the reflectors in combination with the LED sources enables a homogenous and shadow-free illumination field.

Essential features:

- pleasant, shadow less light
- easily movable, keeps the adjusted position safely
- removable handles, sterilizable
- dimmer
- noiseless fan
- adjustable color temperature

Control elements:

1	Handle
2	Regulator for brightness
3	Regulator for Colour temperature
4	Sensor

Startup:

Activate the operating lamp via one of the control panels as described on page 23. The operating light has a proximity sensor (4) on the front that can switch on and off the lamp without touching it.

Light intensity:

The required intensity can be adjusted individually by turning the dimmer (2).

Reduce light intensity:

Turn regulator (2) to the *left* (rear view of the operating lamp) until it reaches the desired intensity.

Increase the light intensity:

Turn regulator (2) to the *right* (rear view of the operating lamp) until it reaches the desired intensity.

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Adjustable range: 4.200 K - 6.000 K

Increase the colour temperature:

→ rotate the regulator (3) anti-clockwise

Reduce the colour temperature:

→ rotate the regulator (3) clockwise

Light field

Size: 7 cm x 14 cm at a distance of 70 cm from the light source and with well defined outlines, in order to avoid risk of patient dazzling.

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Each time the lamp is turned on, the light intensity will be at the level set when it was turned off the time before.

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Control elements:

- 1** Lever
- 2** Handle
- 3** Parabolic reflector
- 4** Lens
- 5** Proximity sensor (on the back side)

The operating light is mounted to the treatment unit on a flexarm.

It is equipped with two LEDs.

Power supply: 17V-24V AC, 50/60Hz,
 22V-35V DC

Rated input: 9VA (maximum)

Essential features:

- pleasant, shadow less light
- easily movable, keeps the adjusted position safely
- removable handles
- Electronic light regulation from maximum to minimum through a microprocessor
- no fan necessary
- sterilizable and easy to remove handles

Startup:

Switch on/off the operating light when dental unit is on and the operating lamp is activated (see page 23):

- by tilting the lever (**1**)

or

- with Proximity sensor (**5**)

To turn the lamp on and off, place your hand close to the sensor, within a maximum distance of 3 cm. When the command is given, an acoustic signal will be heard (1 beep).

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Illumination level:

Each time the lamp is turned on, the light intensity will be at the level memorised when it was turned off the time before.

Changing the light intensity via lever

- To reduce the light intensity, hold the lever (1) pressed to the left/right side until it reaches the desired level of intensity.
- When the lowest level of intensity is reached, an acoustic signal will be heard (1 beep).

Changing the light intensity via proximity sensor

- To reduce the light intensity, remain with your hand near the sensor (5) until the desired level of intensity is reached.
- When the lowest level of intensity is reached, an acoustic signal will be heard (1 beep).
- continue to hold the hand near the sensor (5) to increase the intensity again

Each time the lamp is turned on, the light intensity will be at the level memorised when it was turned off the time before.

Do block the light beam for patients at risk (i.e. with eye conditions). Use correct protection and precautions.

WARNING - reversible glare due to wrong use of LED operating lamp -
A wrong use of the LED operating lamp can lead to temporary visible disturbance, therefore, when turning on or moving of lamp, never direct light field to eyes of patients, yourself or others. Always ensure a distance of approx. 70 cm between lamp and head of patients.

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
Conformity with the hygienic standards as well as good performance of the equipment are warranted by observance with the instructions to the recommendations for disinfection, cleaning and maintenance.

Drilling instruments:

- Every day before the first treatment, let instruments run for approx. 120 seconds to purge all spray tubes.
- Clean and maintain turbine, handpieces and contra-angle attachments.

	IMPORTANT: Please read the manufacturer's operating manual for the drilling instruments very carefully.
---	--

	IMPORTANT: Instrument hoses may not be cleaned in the thermodesinfektor. Temperature load for the tubes is max. 40°C
---	---

Suction handpieces:

Detach handpiece from hose. Clean accessible hose ends and handpieces. Lubricate O-rings with vaseline and reassemble handpiece to the hose. Unsuitable and not permitted are cleaning agents based on quaternary ammonium compositions, cleaning agents containing phenol and aldehyde as well as foam generating agents.

	IMPORTANT: For dosage directions, please refer to the recommendations given for the individual cleaning and disinfection agents.
---	---

Small suction

Small suction should be equipped with disposable suction nozzle.

- Throw away the suction nozzle after treating each patient and replace with a new one before the treatment of next patient.

Big suction

- sterilize Big suction nozzle by high-temperature steam sterilization (134°C, 15 min)

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Cleaning the filter (suction)

Saliva suction

Cleaning suction hose and suction filter

- take out the suction handle from instrument holder
- pull out the suction hose from the suction filter
- take out the filter screen from the filter and rinse with clean water
- install suction handpiece
- install filter screen into the filter

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Cleaning suction system:

- empty a glass of clean water after each patient with small and big suction hose
- Pump/rinse nozzles/suction handpiece with clear water or distilled water for approx. 30 s 1-2 times a day.

Water filter

The unit has water filter to ensure its normal use. After serving a certain time, the dirt will block the filter core and affect the water flow, therefore, it is necessary to clean and replace the filter core.

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
Spittoon bowl

- Remove mushroom cap and filter screen out from the spittoon bowl.
- Clean the mushroom, filter screen, spittoon bowl and spittoon bowl holder with water and disinfect it.
- install the spittoon bowl
- install mushroom and filter screen into the spittoon bowl

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


Turbine exhaust air oil container

- Empty turbine exhaust air oil container (1) and replace cotton wool. The container is located at the bottom side of the dentist unit

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Cleaning and care instructions

What	How	Wherewith	Notes
	<p>● Clean/disinfect upholstery Wipe upholstery, wiping disinfection</p> <ul style="list-style-type: none"> • Headrest upholstery • Backrest upholstery • Seat upholstery 	<p><u>Cleaning:</u></p> <ul style="list-style-type: none"> • Soft sponge, lukewarm water, soap. • Clean heavily soiled areas with a soft brush. <p>After being cleaned, the upholstery should be wiped dry with a cloth.</p> <p><u>Disinfection:</u> We recommend applying alcohol-free disinfectant spray foam for cleaning and disinfecting imitation leather upholstery.</p>	<p>Remove soap residue with a damp sponge and dry it with a soft, lint-free cloth.</p> <p>A cloth is not suitable for cleaning! Never use cleaning agents or household cleaners.</p> <p><u>Discolouration:</u> In most cases, discolouration stains that are found immediately can be removed with the right cleaning tools and cleaner for imitation leather.</p> <p>Discolouration stains that are not detected immediately cannot be removed.</p> <p>Protect the upholstery from direct sunlight.</p>
	<p>● Clean/disinfect upholstery Wipe upholstery, wiping disinfection</p> <ul style="list-style-type: none"> • Backrest upholstery • Seat upholstery 		

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**Keep liquids away from electrical
connections and components!**

- Most electrical control switches ARE NOT hermetically sealed.
- When a switch get wet, it can self-activate, or perhaps short out tad damage other system components!
- Keep switches dry.
- Avoid costly equipment damage or failure.
- Please do not spray disinfectants directly at or into electrical switches and controls.
- Simply spray these solutions onto a cloth (not saturate the cloth) and wipe down the control surfaces.

Any not recommended maintenance/change by the manufacture may cause risk of injury or shock!

- Don't let more than one patient sit in the chair at a time
- Don't stand anywhere on the dental chair
- Don't operate the dental chair around any flammable liquids.

Handpiece

Handpiece is a delicate device, so it shall not be pounded or dropped.

It shall be cleaned and lubricated with special detergent and lubricant every day. As for the specific method, refer to the instruction manual of the respective manufacturer.

Oil spray lubrication shall be conducted before high temperature sterilization.

Disinfection

Cannula of Syringe

- sterilization in steam sterilizer: at 134°C for 3 min

Suction handpiece

- sterilization in steam sterilizer: at 134°C for 3 min

Silicone Mat for instrument tray

- sterilization in steam sterilizer: at 134°C for 3 min

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Overall failure

Failure	Troubleshooting
Failure 1: After power-on, the treatment unit isn't responding	<ol style="list-style-type: none"> 1. Check whether AC 230 power supply input is in good condition. 2. Check whether the 10A fuse at power switch is broken.
Failure 2: After power-on, no water or air is supplied to the treatment unit	<p>Check water source change-over, pressurizing/pressure regulators are in good position;</p> <ul style="list-style-type: none"> • when water source switch is in tap water position, the pressurizing/pressure relief switch should be in OFF-position („Bottle-OFF“); • when water source shift is in bottle water position, the pressurizing/pressure relief switch should be in ON-position („Bottle-ON“).

Common Failures

Failure	Troubleshooting
Failure 1: Operating chair fails to move	<p>Check whether there is 24V AC voltage input at J6 port of chair control board.</p> <ol style="list-style-type: none"> 1. If there is 24V AC voltage input at J6 port, check whether the two yellow keys can operate the movement of the chair; if not, the control board is likely damaged, replace the control board to resolve the problem. 2. If there is no 24V AC voltage input at J6 port, check whether the chair power supply plug on flat cable is in good connection; if not, tighten the screw.
Failure 2: Operating chair moves automatically	<p>In case that operating chair moves automatically, perform troubleshooting from the following aspects:</p> <ol style="list-style-type: none"> 1. Assistant Section Touch Control Board: The assistant section control board may be corroded due to damp absorption or entry of water and thereby causes malfunction of the keys. Inspection method: Disconnect the communication cable terminal on the control board, operate the machine by using pedal and dentist element to check whether the failure remains; if the failure is resolved: it indicates that the control board is damaged and it needs to be replaced; if the failure remains, re-connect the communication cable terminal, and continue with the next step. 2. Dentist element touch control board: The dentist element control board may be corroded due to damp absorption or entry of water and thereby causes malfunction of the keys. Inspection method: Disconnect the communication cable terminal on the control board, operate the machine by using pedal and assistant element to check whether the failure remains; if the failure is resolved, it indicates that the control board is damaged, and it needs to be replaced resolve the problem.
Failure 3: Water heater fails	<ol style="list-style-type: none"> 1. If water heater fails, check whether there is 23V AC voltage output at J13; if yes, the heater may be broken and it needs to be replaced; if no, the control board or the heater is could be damaged; if heater cannot be disconnected, and J14 is always closed, the heater may be broken and it needs to be replaced. 2. If there is no voltage output at J13, check whether fuse F4 is broken; if yes, replace with a 6A fuse.
Failure 4: Handpieces operate normally but no water or air is supplied	<ol style="list-style-type: none"> 1. First check whether the handpiece operates currently in dry grinding mode. For more details, see the operating instructions. 2. Check whether the water/air hose line is in good connection in dentist element. If the hose line is in good condition, then check whether there is 24V DC output at water/air circuit solenoid valves; if yes, the solenoid valve may be broken, and it needs to be replaced; if not, measure voltage output at reserved port J33, if there is a voltage output, replace dentist element control board, if not, inspect the power supply in water unit. 3. foot-operated micro switch: Check whether the foot-operated micro switch resets normally; if not, replace the micro switch. 4. Replace water unit control board.

Failure	Troubleshooting
Failure 5: Program chair fails to store positions	Generally, failure of program chair storage function is a result of a damaged angle sensor, replace the corresponding angle sensor.
Failure 6: There is water or oily dirt in air hose.	Check air supply source and air filter, clean or replace relevant components.

Ultimate Comfort

C20-0087

2302


Anhang


Appendix


Ritter® Concept GmbH
Freiburger Str. 45, D-88400 Biberach

Fon: +49 (0) 7351 52925-10
Fax: +49 (0) 7351 52925-11

e-mail: info@ritterconcept.com
Internet: <http://www.ritterconcept.com>

Instructions about Hydro-group PCB
 1. There are 3 pins at J90 and the reserved connector is placed on the lowest pin. If E. Micromotor is installed, please short circuit the upper 2 pins with the connector.
 2. Default short circuit for OP Light: 2-5 3-6 4-1
 J11 output: with LED-light connection: 0V @Pin A, 12V @Pin C; with Halogen-light connection: 0V @ Pin A, 10V @ Pin B, 12V @Pin C.
 3. LED Indicator status: D4(indication for 24V) and D24(indication for 3.3V): constant ON; D1(indication for CPU running): Blinking; D23 (indication of good COMM with Dentist Table): Blinking.

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Inspection protocol for safety inspection according DIN EN 62353			
Client (dental practice)		Contractor	
Device		Used measurement device	
Device type:	Dentale Behandlungseinheit	Device type:	Electrical Safety Tester
Manufacturer:	Ritter Concept GmbH	Manufacturer:	
Type:		Type:	
Serial No:			
Netzspannung:V	Insulation class:	1
Safety inspection according:	EN 62353/VDE 0751	Type of applied parts:	Typ BF
Visual inspection			
<input type="checkbox"/> Protective earth conductor not damaged (for insulation class 1 only)			
<input type="checkbox"/> Cover and mechanical parts not damaged			
<input type="checkbox"/> Insulating parts not damaged			
<input type="checkbox"/> Device connecting cables including plugs without defects			
<input type="checkbox"/> Labels exist or completed, nameplate correctly			
<input type="checkbox"/> Others			
Measurement			
		measured value	limit value
Protective earth resistance			<0,3Ω
Insulating resistance			>2 MΩ
Test voltage			~500 V
Device leakage current (Equivalent device leakage current)			<10 mA
Encl. leakage current (Housing leakage current)			<5 mA
Leakage current of application part (Equivalent patient leakage current)			<5 mA
Test voltage leakage current (nominal Voltage ±10%)			min.: -10% U _{Netz} max. + 10% U _{Netz}
Functional testing			
<input type="checkbox"/> Function and safety testing without defects			
<input type="checkbox"/> The unit has safety deficiencies:			
<input type="checkbox"/> Short-circuit L-N			
<input type="checkbox"/> Body circuit L/N-SL			
Test cycle: months		
Next test date:		
Tester:			
.....			
Place, Inspection date	Signature		

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Guidance and manufacturers declaration - electromagnetic emission		
The Ultimate Comfort is intended for use in the electromagnetic environment specified below. The customer or the user of the model Ultimate Comfort should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions according CISPR 11	Group 1	The Ultimate Comfort uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The Ultimate Comfort system is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions according CISPR 11	Class B	
Harmonic emissions according IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions according IEC 61000-3-3	Compliant	

Guidance and manufacturers declaration - electromagnetic immunity			
The Ultimate Comfort is intended for use in the electromagnetic environment specified below. The customer or the user of the Ultimate Comfort should assure that it is used in such an environment.			
Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) according IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst according IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge according IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Interruptions and voltage variations on power supply input lines according to IEC 61000-4-11	<5% U_T (>95% dip in U_T for 0.5 cycle) 40% U_T (60% dip in U_T for 5 cycles) 70% U_T (30% dip in U_T for 25 cycles) <5% U_T (>95% dip in U_T for 5 s)	<5% U_T (>95% dip in U_T for 0.5 cycle) 40% U_T (60% dip in U_T for 5 cycles) 70% U_T (30% dip in U_T for 25 cycles) <5% U_T (>95% dip in U_T for 5 s)	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Ultimate Comfort requires continued operation during power mains interruptions, it is recommended that the WED-3000 B Ultrasound Scanner be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field according IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
REMARK U_T is the a. c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity			
The Ultimate Comfort is intended for use in the electromagnetic environment specified below. The customer or the user of the Ultimate Comfort should assure that it is used in such an environment.			
Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the Ultimate Comfort, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance
Conducted RF according IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	$d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$
Radiated RF according IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = \left[\frac{3,5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m) ^b . Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
REMARK 1 At 80 MHz and 800 MHz, the higher frequency range applies. REMARK 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Ultimate Comfort is used exceeds the applicable RF compliance level above, the Ultimate Comfort should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Ultimate Comfort.			
b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the Ultimate Comfort			
The Ultimate Comfort is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Ultimate Comfort can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Ultimate Comfort as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter P in W	Separation distance according to frequency of transmitter l in m		
	150 KHz to 80 MHz $d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.20	1.20	2.30
10	3.80	3.80	7.30
100	12.00	12.00	23.00
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
REMARK 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. REMARK 2: These guidelines may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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Hersteller: ALPRO MEDICAL GmbH

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RiTren

For optimal hygiene of process water in dental unit water lines.

- **RiTren** was developed for use for continuous disinfection of process water in dental unit water lines (incl. instruments and cup filler) for all dental treatment units and it prevents algae formation.
- **RiTren** effective against a broad spectrum of germs
- **RiTren** prevents the formation of lime deposits and (in heating the water and longer service life) - at a water hardness of 2°dH (0.36 mmol/l) up to a water hardness of 20°dH (3.57 mmol/l).
- **RiTren** is due to the greatest purity raw materials of the highest quality, physiologically tested and classified as harmless; It conforms to the requirements of the DAB and the German Cosmetics Regulation.
- **RiTren** is non-corrosive.
- **RiTren** has a very pleasant taste (mint) and promotes the patient's well-being.
- Precondition for the application of RiTren is the Aqua-Bottle-System with 1% RiTren process water.
- **RiTren** has been evaluated for its microbiocidal efficacy by the Institute Fresenius, Stockach 10/95 and 08/96 (RiTren).

Packing unit: 1 l **RiTren**
 Order/Part No: U000003911 capacity: 1 l

We look forward to your order:
 Ritter@Concept GmbH
 Freiburger str. 45 · 88400 · Biberach
 Germany

Fon: +49 7351 52 925-37
 Fax: +49 7351 52 925-21
 info@ritterconcept.com
 www.ritterconcept.com

Manufacturer: ALPRO MEDICAL GmbH

Technical modifications reserved!

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Terms of Guarantee

You are now owner of a high quality Ritter product. Thank you very much for your trust in Ritter Concept GmbH and would like to maintain this! You have the warranty right of 12 months. We point to the warranty terms in our General Terms and Conditions.

We warrant, that our products are free of faults by the time of handover to the supplier. In particular, we warrant, that the production is taking place based on the acknowledged technical rules and the relevant regulations. We also warrant, that the product characteristics are available, which are determined in the contractual documents. In line with the warranty we perform a remedy of defects, when substances, parts or performances are identified as flawed or contrary to contract. The warranty starts with the date of delivery. Deficiencies are to be communicated with Ritter Concept GmbH immediately after receipt of goods and immediately after discovery. The customer is permitted to send back products without prior approval of Ritter Concept GmbH. The customer has to select the most economic way of transport.

No warranty claims exist in the following cases:

- using the product outside of or contrary to its intended use, or operating the product outside of its specifications
- failure to observe the instructions on use, maintenance or installation
- improper handling or neglected maintenance
- use of non-approved spare parts or accessories
- product modifications or conversions performed by customer or third parties
- usual wear and tear or other changes to or deterioration of the product attributable to:
 - external influences (e. g. knocks, blows, mechanical impacts)
 - water, fire
 - improper method of storage, handling or installation
 - unusual ambient conditions, especially conditions upon receiving or operating conditions at the installation location or force majeure
 - calcification or corrosion resulting from neglected maintenance or improper handling
 - contaminants in the air and water supply
 - unusual or impermissible chemical or electrical influences
 - soiling, discolouration or wear and tear on padding

The warranty as a rule excludes parts subject to normal wear, especially lamps, glass and rubber components and cables as well as plastic parts and surfaces regarding their colour-fastness.

All of our obligations under the warranty will be extinguished if the delivered goods are modified, improperly handled or processed. We do not warrant co-supplied or integrated third-party products. However, if requested we will assign to the customer our warranty claims against the vendor of the third-party product.

In order to activate the warranty, it is necessary to send a copy of the record of delivery and the warranty card to us following installation.