# Leadex 70 DC

Intraoraler Röntgenstrahler

intraoral x-ray system at constant potential

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
This operating manual should always be kept handy near the installation.

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THE RADIOGRAPHIC SYSTEM DESCRIBED IN THIS MANUAL REFERS TO A WALL INSTALLATION.

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THE MANUFACTURER RESERVES ITSELF THE RIGHT TO MODIFY THE PRODUCTION AND THE MANUAL WITHOUT NOTICE. COPIES, EVEN PARTIAL, OF THIS MANUAL ARE ALLOWED ONLY FOR IN-HOUSE USE.

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THE MANUFACTURER SHALL NOT BE LIABLE FOR AN INCORRECT USE OF THE INFORMATION CONTAINED IN THIS MANUAL.

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**PRELIMINARY INFORMATION**

Before beginning to use the “**Leadex70 DC**” radiographic system, it is advisable to carefully read and follow the instructions contained herein, so as to obtain the best possible performance.

Always pay close attention to the

**CAUTION  
WARNING  
PLEASE NOTE**

messages when operating the system.

**LEGEND**

**CAUTION**

*The word **CAUTION** identifies those occurrences which might compromise the operator's personal safety or cause injuries to people.*

**WARNING**

*The word **WARNING** identifies those occurrences which might compromise the radiographic system's performance.*

**PLEASE NOTE**

**PLEASE NOTE** serve to give special indications so as to facilitate maintenance or make important information clearer.

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**INFORMATION FOR THE USER**

Dear Customer,  
thanks for having chosen the „**Leadex70 DC**“ radiographic system.  
It is designed and manufactured by „de Götzen® S.r.l.“ and is the result of many years of experience in the field of radiology and in the application of advanced electronics.

This high performing system represents a further development of technological research at the service of dental radiography.

**PLEASE NOTE**

This manual does not contain all the recommendations and the obligations relative to the possession of a source of ionising radiations – since they do vary from State to State – but only the most common ones.  
The user must consult its country's legislation so as to fulfil all local obligations.

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Inappropriate use or any arbitrary tampering with, exempt the manufacturer of the „Leadex70 DC“ radiographic system, from any service under warranty or from any other liability.

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The warranty is valid only if the following precautions are taken:

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- any repairs, modifications, adjustments, recalibrations must be performed only by the manufacturer
- the installation must be made by professionally qualified technicians according to the regulations in force
- the system must be installed and used in compliance with the instructions given in this Manual and for the purposes and applications for which it was designed
- the power supply must be adequate to supply the required power indicated in the radiographic system’s nameplate data
- in order to safeguard one’s warranty rights, please fill in the enclosed Warranty Document, immediately after the installation is completed, together with the technician

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All claims for damage or miscarriage regarding the shipment must be pointed out in the presence of the shipping agent.

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In case of miscarriages, or actual or suspected damage, the receiver shall indicate the proper reserves on the way-bill or on the consignment note.

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A few safety recommendations are listed here below which should be followed when using the „Leadex70 DC“ radiographic system.

**CAUTION****PROTECTION AGAINST RADIATIONS**

*"The general principles regarding safety and protection of workers and people" must always be applied when using the unit:*

1. *Justification of the practice*
2. *Protection Optimisation (ALARA)*
3. *Reduction of the limits of individual dose and risks*

*The radiographic system must only be used by authorized and qualified personnel.*

*All personnel present during the radiographic examination must comply with safety measures foreseen in regards to protection from radiations. For his own safety, the operator must always keep a distance of more than 2 meters from the radiographic unit.*

*To protect the patient from unnecessary exposure to radiations, additional anti-radiation protections may be used whenever necessary (aprons, collars, etc...)*

**ATTENTION X-RAYS****CAUTION****MECHANICAL RISK**

*Before removing the tubehead from the positioning arm, RELEASE THE SPRING because the joint might burst open and hit the operator.*

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**CAUTION****ELECTRIC SAFETY**

*The radiographic system contains high voltage.*

*When inspecting internal parts, always turn off the power before touching any electric part.*

*The unit must be used only in environments that are in compliance with all electric safety standards set forth for medical environments.*

*The unit is NOT equipped with protections against penetration of liquids; it will therefore be necessary to make sure that no water or other liquids penetrate inside so as to avoid short circuits or corrosion.*

*Always disconnect the radiographic system from the power supply before beginning cleaning and disinfecting operations.*

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**CAUTION****PROTECTION AGAINST EXPLOSIONS**

*The radiographic system MUST NOT be used in the presence of disinfectant, flammable or potentially explosive gases or vapours that might catch fire and cause damage.*

*In case these disinfectants have to be used let the vapour completely disperse before turning on the radiographic system.*

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- ▶ Radiation Performance Standard 21 CFR, Subchapter J

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**RADIOGRAPHIC SYSTEM**

The radiographic system „**Leadex70 DC**“ guarantees the maximum safety both for the operator and the patient.

It is built in compliance with the following European Directives:

- ▶ 93/42/EEC  
MEDICAL DEVICES
- ▶ 73/23/EEC  
LOW VOLTAGE
- ▶ 89/336/EEC  
ELECTROMAGNETIC COMPATIBILITY
- ▶ EURATOM 96/29  
IONISING RADIATIONS

and in compliance with the following American Standard:

- ▶ American Radiation Performance Standard 21 CFR, Subchapter J

The following protective measures were adopted in the design and construction of the unit:

- protection against the risk of electric injuries, ensured by a grounded protection conductor
- protection against leakage radiation, made negligible by the shielded casing
- protection against excessive radiations, thanks to the immediate activation of the safety device
- protection against continuous service, since the system is designed, according to standards, so as not to allow use in radioscopy
- protection for the patient against dangerous radiations, obtained by means of the high frequency technology capable of producing a constant and hard radiation

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"ELECTRO-MEDICAL" CLASSIFICATION

According to paragraph §5 of the general safety regulations EC EN 60 601-1/1998 on safety of medical equipment, the system is classified as: **Class I - Type B**

"MEDICAL DEVICES" CLASSIFICATION

According to the classification rules indicated in attachment IX of the EEC Directive 93/42 on medical devices the system is classified as: **Class IIb**

"E.M.C." CLASSIFICATION

According to paragraph §4 of the EEC EN 55011, the system is classified as: **Group 1 – Class B**

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The „Leadex70 DC“ radiographic system (Fig. 1) consists of:

### 1. TIMER

The timer is the control panel used to manage the times and to safely use the tubehead.

For making the exposure it is available the control button with safety button.

The timer can be connected to n° 2 tubeheads.

The timer is „multi-technological“ type and it is able to control tubehead of different technologies: alternate current tubehead and direct current tubeheads.

In case of alternate current tubeheads the technology of the timer is „self-compensating“.

Depending on the line voltage fluctuation the microprocessor automatically modifies the predetermined exposure time guaranteeing a constant dose to the patient.

This technological expedient avoids the repetition of the exposure because of over/under exposure errors.

### 2. PANTOGRAPH TYPE ARM

Thanks to the new shape and new mechanisms of the positioning arm, it can be adjusted in height and depth so as to precisely explore any spot in its reach.

It is made of light alloy with an ABS coating.

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### 3. TUBEHEAD

The intra-oral tubehead is of the monoblock type and the light alloy housing is divided into two compartments.

The high voltage transformer, the x-ray tube and the expansion chamber are submerged in highly dielectric.

The expansion chamber guarantees an adequate compensation to oil expansion for the entire temperature range.

The x-ray tube is located in the back part of the container this allowing a source-skin distance that is 50% higher than traditional structures.

In the second compartment are placed the main electronic board and the control electronic board.

### 4. CONE

Made of the transparent polycarbonate, it allows for:

- correct focal spot to skin distance
- dimension, direction and centering of x-ray beam
- realization of different radiographic technique (bisecting and parallel technique)

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

## NAMEPLATES

The identification tags on the tubehead, on the timer and on the cone indicate the model number, the serial number, the manufacturing date and the symbols of the main technical characteristics.

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


**CONFIGURATION**

The radiographic system is factory configured "standard mode".

*On the control panel the led relevant to the following exposure parameters will light up:*

Nummer des gewählten Röntgenkopfes LED 1		No of the selected tubehead LED 1
Gelieferter Tubus LED 8" = KURZTUBUS LED 12" = LANGTUBUS		supplied cone LED 8" = SHORT CONE LED 12" = LONG CONE
Art des Röntgenkopfes LED AC = WECHSELSTROM LED DC = GLEICHSTROM		type tubehead LED AC = ALTERNATE CURRENT LED DC = DIREKCT CURRENT
Röhrenspannung LED 70kV		radiographic voltage LED 70 kV
Röhrenstrom LED 8mA		radiographic current LED 8 mA
Patiententyp LED ERWACHSENER		type of patient LED ADULT
Aufnahmetechnik KONVENTIONELL LED D		radiographic technique CONVENTIONAL LED D

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*The following exposure times (sec) have been stored:*

**0,020 – 0,025 – 0,032 – 0,040 – 0,050 – 0,063 – 0,080 – 0,100 – 0,125 – 0,160  
– 0,200 – 0,250 – 0,320 – 0,400 – 0,500 – 0,630 – 0,800 – 1,00 – 1,250 – 1,600  
– 2,000 – 2,500 – 3,200**

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**PLEASE NOTE**

These times are in compliance with current CEI EN 60601-2-7 (1999) norms and with the ISO 497 series R'10 recommendations.

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*Certain exposure values have been predefined which depend on the selection of the operating parameters:*

- ▶ cone (8"/12")
- ▶ type of patient (ADULT/CHILD)
- ▶ radiographic technique
- ▶ intra-oral test

**PLEASE NOTE**

If one so desires, it is possible to change this values. (refer to Charter 5 and 6)

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*Possible modifications of the exposure values*

- ▶ radiographic voltage (60 kV/70 kV)
- ▶ radiographic current (4 mA/8mA)
- ▶ type of patient (ADULT/CHILD)
- ▶ radiographic technique

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*Possible modifications of the exposure values*


- ▶ cone (8"/12")
- ▶ No of tubehead
- ▶ No control button

inside the timer, by changing the dip-switch position

THIS OPERATION MUST BE CARRIED OUT BY THE INSTALLER ONLY!

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Bring the main switch located on the upper part of the timer to the "I" position (ON)

Bring the key switch to the "I" position (ON)

1. the green light turns on indicating that the system is powered
2. the LEDs of the set parameters automatically light up
3. the exposure time is shown on the display

**CAUTION**

**If an error is detected when the system is turned on, the anomaly is indicated as follows:**

- an intermittent sounds off (beep)
- intermittently turn on  **MALFUNCTIONING INDICATOR**
- the error code (E ...) appears on the display (refer to Chapter 8)
- all control panel functions are inhibited

**In this case turn off the timer and then turn it back on.  
If the error should repeat itself, call the "Assistance Service".**

**PLEASE NOTE**

The exposure time and parameters which appear on the display are the last that were set before the timer was turned off. If the timer remains inactive for a few minutes, it switches to the stand-by mode. Press any key on the control panel to restore it to the operative mode.

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## CHECK THE SELECTED PARAMETER - 2°

Before making the exposure check that the control panel selected parameter (from Step 1 to Step 8) are suitable for the radiographic exam.

STEP 1  
check the selected tubehead

The LED of the desired tubehead should be turned on

### LED Rx 1 ON

indicates that is selected the tubehead connected to the timer XRAY1 terminal block

### LED Rx 2 ON

indicates that is selected the tubehead connected to the timer XRAY2 terminal block

*to change the selection press again the button*

STEP 2  
check the selected radiographic distance CONE

The LED of the cone length (source-skin distance = SSD) in use should be turned on

### LED 8" ON

indicates that the selected tubehead is equipped with 8" = 20cm (SSD) cone

### LED 12" ON

indicates that the selected tubehead is equipped with 12" = 31cm (SSD) cone

*to change the selection call the "Assistance Service"*

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The LED of the type of selected tubehead should be turned on

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**LED DC ON**

indicates that the selected tubehead works in direct current technology

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The LED of the desired radiographic voltage should be turned on

**STEP 4**  
check the  
selected  
radiographic  
voltage

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**LED 60 kV ON**

indicate that the radiographic system is set with the high contrast radiodiagnostic technology

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**LED 70 kV ON**

indicates that the radiographic system is set with the low contrast radiodiagnostic technology

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**PLEASE NOTE**

60 kV mode can be selected with the "Leadex70 dc" radiographic system only.

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The LED of the desired radiographic current should be turned on

### STEP 5

check the selected radiographic current

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#### LED 4 mA ON

indicates that the radiographic system is set with reduced dose. It is advisable the use of digital radiographic technic.

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#### LED 8 mA ON

indicates that the radiographic system is set with nominal dose. It is advisable the use of conventional radiographic technique.

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#### PLEASE NOTE

The 4 mA mode can be selected with the "Leadex70 DC" radiographic system only.

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The LED of the desired type of patient should be turned on

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#### LED CHILD ON

indicates that the radiographic system is set for a patient with a small physique

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#### LED ADULT ON

indicates that the radiographic system is set for a patient with a large physique

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*to change the selection press the button again*

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## STEP 7

check the selected radiographic technique

### CONVENTIONAL TECHNIQUE (FILM)

The LED of the desired speed film should be turned on

#### LED **D ON**

radiographic system is set for use with D speed film

#### LED **E ON**

radiographic system is set for use with E speed film

#### LED **F ON**

radiographic system is set for use with F speed film

*to change the selection press the button again for 3 s until the acoustic signal*

#### **PLEASE NOTE**

With films it is advisable to use a radiographic current of 8 mA.  
(refer to STEP 5)

#### **PLEASE NOTE**

After the modification, default exposure values will be automatically changed.

### DIGITAL TECHNIQUE (SENSOR)

The LED should be turned on

*to change the selection press the button again for 3 s until the acoustic signal*

#### **PLEASE NOTE**

With films it is advisable to use a radiographic current of 8 mA.  
(refer to STEP 5)

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### PERIAPICAL EXAM

The LED of the selected teeth should be turned on

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*to change the selection press the key relative to the desired tooth*

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### OCCLUSAL EXAM

The LED of the selected type of test should be turned on

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#### **LED MANDIBULA ON**

radiographic system is set for occlusal exam of the lower jaw

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#### **LED MAXILLA ON**

radiographic system is set for occlusal exam of the upper jaw

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#### **LED ANT ON**

radiographic system is set for anterior bite-wing exam

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#### **LED MAXILLA ON**

radiographic system is set for posterior bite-wing exam

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*to change the selection press the key relative to the desired exam*

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### POSITIONING THE PATIENT - 3°

Following the standard intra-oral procedures positioning the patient.

### 4° - POSITIONING FILM or SENSOR

Positioning either the film or the digital sensor depending on the technique to be used.

### 5° - POSITIONING CONE

Following the standard intra-oral procedures positioning the cone.

### 6° - CHECK ON THE DISPLAY THE SELECTED TIME

Before proceeding with the exposure check on the display the selected time

*to change the selection press the following keys*

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#### **WARNUNG**

*This modification brought to the exposure time is momentary and will be lost unless it is saved (refer to Chapter 6).*

*To restore the previous values, press one of the keys with the led turned off on the control panel.*

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## 7° - MAKE THE EXPOSURE

1. Take the control button of the timer relevant to the selected tubehead and keep a safety distance of a least (2 meters) from the tubehead, so as to be able to constantly check the radiographic exposure
2. Advise the patient to remain still
3. On the control button press the  X-RAY key and keep it pressed until the acoustic signal (beep) stops and led yellow turn off 

### PLEASE NOTE

If the "X-RAY" key is released early, the exposure is immediately interrupted and the E12 error message appears on the display.

4. At the end of the exposure the green led intermittently turn on  PAUSE
5. The display indicates the actual exposure time
6. All the timer functions are inhibited

### PLEASE NOTE

The pause time is necessary to allow the x-ray tube to cool down.

This time is calculated by the microprocessor, depending on the exposure time, with a ratio of 1:32 (32 second of pause are required for each second of exposure)

### **A NEW EXPOSURE WILL BE POSSIBLE**

**AFTER THE GREEN LED  HAS TURNED OFF**

**REPEAT THE OPERATIVE SEQUENCE FROM POINT 2 TO POINT 7**

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**CHART OF DEFAULT  
EXPOSURE VALUES**

The chart indicates the "Leadex70 DC" radiographic system predefined exposure values.

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
**PLEASE NOTE**

To modify the default exposure times. (*refer to Chapter 6*)


**LANGTUBUS12"  
(SSD = 31 cm)**

**KONVENTIONELLE RÖNTGENTECHNIK (FILM)**

**12" CONE LONG  
(SSD = 31 cm)**

**CONVENTIONAL RADIOGRAPHIC TECHNIQUE  
(FILM)**

ERWACHSENER/ADULT																								
VORPROGRAMMIERTE BELICHTUNGSZEITEN (s)	0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200	
PROGRAMMED EXPOSURE TIMES (s)	0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200	
<b>70 kV - 8 mA</b>																								
FILM D												I		M	Op									
FILM E											I		M	Op										
FILM F									I		M		Op											

ERWACHSENER/ADULT																								
VORPROGRAMMIERTE BELICHTUNGSZEITEN (s)	0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200	
PROGRAMMED EXPOSURE TIMES (sec)	0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200	
<b>60 kV - 8 mA</b>																								
FILM D															I		M	Op						
FILM E												I		M	Op									
FILM F											I		M	Op										


KIND/CHILD																								
VORPROGRAMMIERTE BELICHTUNGSZEITEN (s)	0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200	
PROGRAMMED EXPOSURE TIMES (sec)	0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200	
<b>70 kV - 8 mA</b>																								
FILM D												I		M	Op									
FILM E											I		M	Op										
FILM F									I		M	Op												


KIND/CHILD																								
VORPROGRAMMIERTE BELICHTUNGSZEITEN (s)	0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200	
PROGRAMMED EXPOSURE TIMES (sec)	0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200	
<b>60 kV - 8 mA</b>																								
FILM D															I		M	Op						
FILM E												I		M	Op									
FILM F											I		M	Op										

## LANGTUBUS12" (SSD = 31 cm)

## 12" CONE LONG (SSD = 31 cm)

**DIGITALE RÖNTGENTECHNIK  
(SENSOR)**

**DIGITAL RADIOGRAPHIC TECHNIQUE  
(SENSOR)**

ERWACHSENER/ADULT																								
VORPROGRAMMIERTE BELICHTUNGSZEITEN (s)	0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200	
PROGRAMMIERTE EXPOSURE TIMES (s)																								
<b>70 kV - 4 mA</b>																								
MAXILLA						I		M	Op															
MANDIBLE					I		M	-	Os															

ERWACHSENER/ADULT																								
VORPROGRAMMIERTE BELICHTUNGSZEITEN (s)	0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200	
PROGRAMMIERTE EXPOSURE TIMES (s)																								
<b>60 kV - 4 mA</b>																								
MAXILLA									I		M	Op												
MANDIBLE								I		M	-	Os												


KIND/CHILD																								
VORPROGRAMMIERTE BELICHTUNGSZEITEN (s)	0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200	
PROGRAMMIERTE EXPOSURE TIMES (s)																								
<b>70 kV - 4 mA</b>																								
MAXILLA					I		M	Op																
MANDIBLE				I		M	-	Os																

KIND/CHILD																								
VORPROGRAMMIERTE BELICHTUNGSZEITEN (s)	0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200	
PROGRAMMIERTE EXPOSURE TIMES (s)																								
<b>60 kV - 4 mA</b>																								
MAXILLA								I		M	Op													
MANDIBLE							I		M	-	Os													


## KURZTUBUS 8" (SSD = 20 cm)

KONVENTIONELLE RÖNTGENTECHNIK (FILM)

## 8" SHORT CONE (SSD = 20 cm)

CONVENTIONAL RADIOGRAPHIC TECHNIQUE  
(FILM)


		ERWACHSENER/ADULT																								
VORPROGRAMMIERTE BELICHTUNGSZEITEN (s)		0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200		
PROGRAMMIERT EXPOSURE TIMES (s)																										
70 kV - 8 mA																										
FILM D										I		M	Op													
FILM E								I		M	Op															
FILM F						I		M	-	Oa																

		ERWACHSENER/ADULT																								
VORPROGRAMMIERTE BELICHTUNGSZEITEN (s)		0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200		
PROGRAMMIERT EXPOSURE TIMES (s)																										
60 kV - 8 mA																										
FILM D												I		M	Op											
FILM E										I		M	Op													
FILM F									I		M	-	Oa													

		KIND/CHILD																								
VORPROGRAMMIERTE BELICHTUNGSZEITEN (s)		0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200		
PROGRAMMIERT EXPOSURE TIMES (s)																										
70 kV - 8 mA																										
FILM D									I		M	Op														
FILM E							I		M	Op																
FILM F					I		M	-	Oa																	

		KIND/CHILD																								
VORPROGRAMMIERTE BELICHTUNGSZEITEN (s)		0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200		
PROGRAMMIERT EXPOSURE TIMES (s)																										
60 kV - 8 mA																										
FILM D												I		M	Op											
FILM E									I		M	Op														
FILM F								I		M	-	Oa														


## KURZTUBUS 8" (SSD = 20 cm)

DIGITALE RÖNTGENTECHNIK (SENSOR)

## 8" SHORT CONE (SSD = 20 cm)

DIGITAL RADIOGRAPHIC TECHNIQUE (SENSOR)

ERWACHSENER/ADULT																								
VORPROGRAMMIERTE BELICHTUNGSZEITEN (s)	0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200	
PROGRAMMIERTE EXPOSURE TIMES (s)	0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200	
70 kV - 8 mA																								
FILM D									I		M	Op												
FILM E								I		M	Op													
FILM F					I		M	-	Oa															

ERWACHSENER/ADULT																								
VORPROGRAMMIERTE BELICHTUNGSZEITEN (s)	0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200	
PROGRAMMIERTE EXPOSURE TIMES (s)	0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200	
60 kV - 8 mA																								
FILM D												I		M	Op									
FILM E										I		M	Op											
FILM F								I		M	-	Oa												

KIND/CHILD																								
VORPROGRAMMIERTE BELICHTUNGSZEITEN (s)	0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200	
PROGRAMMIERTE EXPOSURE TIMES (s)	0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200	
70 kV - 8 mA																								
FILM D									I		M	Op												
FILM E								I		M	Op													
FILM F				I		M	-	Oa																

KIND/CHILD																								
VORPROGRAMMIERTE BELICHTUNGSZEITEN (s)	0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200	
PROGRAMMIERTE EXPOSURE TIMES (s)	0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200	
60 kV - 8 mA																								
FILM D												I		M	Op									
FILM E										I		M	Op											
FILM F								I		M	-	Oa												

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**! WARNING**  
The 23 programmed exposure times may not be modified in the "Leadex70 DC" radiographic system.  
  
Meanwhile you can customize the default exposure values.  
(refer to Chapter 3)

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**! WARNING**  
After customizing the "Chart of default exposure values" (refer to Chapter 5) are not valid any more.

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*For programming the new exposure values press the following keys*

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**PLEASE NOTE**  
The "repeat" function automatically sets in when the key is kept pressed so the time shown on the display scrolls faster.

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*For confirming the new program check the LED of key*

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**LED MEMO ON**  
indicates that it is possible to save the new default exposure value.

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Press the button for 3sec until the acoustic signal to save the new default exposure values.

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**LED MEMO OFF**  
indicates that it is not possible to save the new default exposure value

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**PLEASE NOTE**  
It is not possible to save data when the "range of exposure field" exceeds the programmed exposure time limits.  
(refer to example at the next page)

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THE RANGE OF EXPOSURE FIELD HAS BEEN  
REDUCED BY TWO STEPS

		ERWACHSEN ER/ADULT																							
PROGRAMMED EXPOSURE TIMES (s)		0,020	0,025	0,032	0,040	0,050	0,063	0,080	0,100	0,125	0,160	0,200	0,250	0,320	0,400	0,500	0,630	0,800	1,000	1,250	1,600	2,000	2,500	3,200	
FILM D												I		M	Op										
FILM D											I		M	-	Oa										
FILM E										I		M	Op												
FILM E										I		M	-	Oa											
FILM F										I		M	Op												
FILM F										I		M	-	Oa											

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ON

#### RESTORING ORIGINAL VALUES

- Turn the timer off
- Keep the key pressed turn the timer on
- On the display appears
- Release the key
- Press again the key
- Appears on the display
- Turn off and on the timer

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With the "Leadex70 DC" radiographic system it is possible to visualise certain functional parameters.

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To visualise them proceed as follows:

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1. press simultaneously and keep pressed the keys

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2. press the key associated to the parameter one wishes to:

ANGEZEIGTER PARAMETER	TASTE KEY	DISPLAY PARAMETER
NENNSPANNUNG RÖNTGENSYSTEM		RADIOGRAPHIC SYSTEM NOMINAL VOLTAGE
NETZSPANNUNG		LINE VOLTAGE
SOFTWARE VERSION		SOFTWARE VERSION

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**ERROR MESSAGES**

The following chart gives a list of error messages that may appear while "Leadex70 DC" radiographic system is working.

The chart also includes the causes of the error messages and what to do to solve them.

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<b>ERROR MESSAGES</b>	<b>CAUSE</b>	<b>SOLUTION</b>
E00	RX1 TUBEHEAD IS NOT CONNECTED OR IS OUT OF ORDER	CALL THE "ASSISTANCE SERVICE"
E01	RX2 TUBEHEAD IS NOT CONNECTED OR IS OUT OF ORDER	CALL THE "ASSISTANCE SERVICE"
E02	CORRUPTED EEPROM DATA	CALL THE "ASSISTANCE SERVICE"
E03	EEPROM DATA NOT SAVED PROPERLY	CALL THE "ASSISTANCE SERVICE"
E05	LINE VOLTAGE VALUE NOT INCLUDED WITHIN THE SET LIMITS	CALL THE "ASSISTANCE SERVICE"
E07	LINE VOLTAGE VALUE NOT INCLUDED WITHIN THE $\pm 15\%$ NOMINAL VALUE	CALL THE "ASSISTANCE SERVICE"
E08	THE X-RAY KEY ALWAYS SEEMS TO BE PRESSED	MAKE SURE IT IS NOT JAMMED
E09	ANOMALY IN THE CONTROL PANEL	CALL THE "ASSISTANCE SERVICE"
E12	THE EXPOSURE HAS BEEN PREMATURELY INTERRUPTED	KEEP THE X-RAY KEY PRESSED TILL THE END OF THE EXPOSURE
E20	ANOMALY IN THE TRIAC/RELAY	CALL THE "ASSISTANCE SERVICE"
E21	ANOMALY IN THE ELECTRONIC CIRCUIT	CALL THE "ASSISTANCE SERVICE"
E22	ANOMALY IN THE CONTROL CIRCUIT	CALL THE "ASSISTANCE SERVICE"
E23	INCORRECT DIP-SWITCH CONFIGURATION	CALL THE "ASSISTANCE SERVICE"
E24	THE CONTROL BUTTON DOES NOT CORRESPOND TO THE SELECTED TUBEHEAD	CALL THE "ASSISTANCE SERVICE"
E30	THE TUBEHEAD DOES NOT WORK PROPERLY	CALL THE "ASSISTANCE SERVICE"
E32	THE TUBEHEAD IS NOT IN THE CORRECT MODE	CALL THE "ASSISTANCE SERVICE"
E33	THE TUBEHEAD HAS NOT COMPLETED THE EXPOSURE	REPEAT THE EXPOSURE CALL THE "ASSISTANCE SERVICE"
E40	PROBLEM IN THE FREQUENCY OR REGULATION	CALL THE "ASSISTANCE SERVICE"
E41	THE TUBEHEAD IS NOT CALIBRATED	CARRY OUT WITH THE CALIBRATION
E42	EEPROM DATA NOT SAVED PROPERLY	CALL THE "ASSISTANCE SERVICE"
E43	CORRUPTED EEPROM DATA	CALL THE "ASSISTANCE SERVICE"
E44	OVERVOLTAGE ERROR	CALL THE "ASSISTANCE SERVICE"
E45	ANODE VOLTAGE OUT OF TOLERANCE	CALL THE "ASSISTANCE SERVICE"
E46	ANODE CURRENT OUT OF TOLERANCE	CALL THE "ASSISTANCE SERVICE"
E47	CONTROL CONNECTOR	CALL THE "ASSISTANCE SERVICE"
E48	PROBLEM IN THE REFERENCE VOLTAGE	CALL THE "ASSISTANCE SERVICE"
ERR	MAJOR ERROR	ALL FUNCTIONS ARE DISABLED CALL THE "ASSISTANCE SERVICE"

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#### TECHICAL DATA OF THE RADIOGRAPHIC SYSTEM

GENERATOR	Gleichspannung / at constant potential	GENERATOR
NOMINALE HOCHSPANNUNG	<b>60 kV 70 kV</b>	NOMINAL HIGH VOLTAGE
NENNSTROM	<b>4 mA 8 mA</b>	NOMINAL CURRENT
NOMINALE ELEKTRISCHE LEISTUNG BEI 0,1SEC	<b>560 W (70 kV, 8 mA)</b> <b>480 W (60 kV, 8 mA)</b> <b>280 W (70 kV, 4 mA)</b> <b>240 W (60 kV, 4 mA)</b>	NOMINAL ELECTRIC POWER AT 0,1sec
VERGLEICHSSSTROM-ZEIT-PRODUKT	<b>0,8m As 8 mA 0,1 s</b> <b>0,4 mAs 4 mA 0,1 s</b>	REFERENCE CURRENT-TIME PRODUCT
STRAHLUNGSINTENSITÄT IN DER LUFT	<b>&gt;38 µGy/mAs in 1 m Entfernung vom Brennfleck</b> <b>at 1 m away from focal spot</b>	INTENSITY OF RADIATION IN THE AIR
GESAMTFILTERUNG	<b>2 mm Al bei/at 70 kV</b>	TOTAL FILTRATION
HALBWERTSSCHICHT BEI 70=KV	<b>&gt;1,6 mm Al</b>	HALF VALUE LAYER (HLV) AT 70kV
STREUSTRABLUNG	<b>&lt;0,25 mGy/h in 1 m Entfernung</b> <b>&lt;0,25 mGy/h to 1 m away from focal spot</b>	LEAKAGE RADIATION
LINEARITÄT	<b>10%</b>	LINEARITY
REPRODUZIERBARKEIT	<b>0,05</b>	RIPRODUCIBILITY
ELEKTRISCHE KLASSIFIZIERUNG	<b>Klasse I - Typ B - intermittierender Betrieb</b> <b>Class I - Type B - intermittent service</b>	ELECTRICAL CLASSIFICATION

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## ACCURACY OF TECHNICAL DATA

NOMINAL VOLTAGE OF X-RAY TUBE
NOMINAL CURRENT OF THE X-RAY TUBE
SELECTED EXPOSURE TIME

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## ENVIRONMENTAL CHARACTERISTICS

OPERATIVE TEMPERATURE
WAREHOUSE TEMPERATURE
HUMIDITY

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
QUELLE-HAUT ABSTAND	
KURZTUBUS	<b>20 cm (8")</b>
LANGTUBUS	<b>31 cm (12")</b>
RECHTECKTUBUS	<b>31 cm (12")</b>
DURCHMESSER DES RÖNTGENSTRAHLS	
KURZTUBUS	<b>≤60 mm</b>
LANGTUBUS	<b>≤60 mm</b>
RECHTECKTUBUS	<b>44x35 mm</b>

## TECHNICAL DATA OF THE CONE

SOURCE-SKIN DISTANCE (SSD)
SHORT CONE
LONG CONE
RECTANGULAR CONE
DISTANCE OF X-RAY BEAM
SHORT CONE
LONG CONE
RECTANGULAR CONE

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
**1,2 kJ/min  
(1,8 kUT/min)**

**ERHITZUNGSKURVE / HEATING CURVE**


## THERMAL CHARACTERISTICS OF THE TUBEHEAD

TUBEHEAD'S HEAT ACCUMULATION CAPACITY

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**(70 kV 8 mA) ±10%**

NOMINAL HIGH VOLTAGE and  
MAXIMUM CURRENT

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**0,8 mm Al bei/at 70 kV**

TUBE INHERENT FILTRATION

ANODENMATERIAL

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ANODE MATEERIAL

ANODENNEIGUNG

**20°**

ANODE INCLINATION

ANODENWÄRMEBELASTUNG

**7 kJ (10 kUT)**

ANODE HEAT LOAD

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**17,5 W**


MAXIMUM CONTINUOUS HEAT DISSIPATION

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In order to guarantee safety of the radiographic system, it is necessary to set up a maintenance schedule.

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The owner is responsible for organising and observing a maintenance schedule which must be executed by qualified technicians who must be able to certify their work with a "Conformity Declaration".

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**CAUTION**

***Run an inspection on the radiographic system and on its operation when it is installed and every twelve months.***

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Once a year, lubricate the pins and bushes of the wall plate and the positioning arm, as specified. (refer to INSTALLATION & MAINTENANCE MANUAL)

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**WARNING**

***Do not lose the adjustment key that comes with the system, since, in time, it could become necessary to make readjustments.***

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**WARNING**

***If the parts should become hard to move or should squeak, call the "Assistance Service".***

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**CLEANING THE OUTER SURFACE**

Use a soft cloth dampened with water and soap to clean the outer surfaces.

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The spacer cone may be cleaned with cotton wool soaked with surgical alcohol.

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In case of a malfunction, send the defective part using the original packaging to:

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### SCRAPPING

The use of the WEEE symbol  indicates that this product may not be treated as household waste, but must be treated separately, in conformity to the Directive 2002/96/CE.

By ensuring this product is disposed of correctly, you will help to protect the environment.

For more detailed information about the recycling of this product, please contact your local authority, your house waste disposal service provider or where you purchased the product.

#### **CAUTION**


**To avoid any risk of environmental contamination, don't dispose the device and its accessories together with the house waste materials.**

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**CHAPTER  
12****ATTACHMENTS**

The manufacturer undertakes to supply, upon request, drawings, circuit diagrams, component parts lists, instructions or other information needed by qualified technical personnel to perform repairs on those parts of the "Leadex70 DC" radiographic system which may be repaired.

# Leadex70 DC

---

distributed by **Ritter® Concept GmbH**

Bahnhofstraße 65  
08297 Zwönitz  
GERMANY

e-mail: [info@ritterconcept.com](mailto:info@ritterconcept.com)

made by **de Götzen® S.r.l.**

Via Roma 45  
21057 OLGiate OLONA VA  
ITALY

e-mail: [degotzen@degotzen.com](mailto:degotzen@degotzen.com)