Contact World Contact Lite Comfort

C20-0020
2202

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Dear user!

Our congratulations on the choice of your new Ritter product and many thanks for your confidence! This modern device was produced by our development- and production engineers with great carefulness.

You and your team will gain absolute joy and benefit of your new Ritter device, if you and your employees, after proper mounting on site, follow the technical and application-oriented advices from the handing over demonstration.

Despite of a careful introduction into the operation, smaller or greater questions will arise in the course of time. This operating manual shall give you answers to these questions. To make day-to-day work easy for you and to be sure, that all technical or legal regulations and the maintenance notes receive their attention, we recommend a careful examination of this manual for you and your responsible employees too.

If you face technical problems, our service hotline
+49 7351 52 925-37 is available for you.
And now a lot of joy with your Ritter equipment.

Abb.: Dentale Behandlungseinheit "Contact World H"
Fig.: Dental treatment unit "Contact World H"

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No.	Description
1	Assistant element
2	Water unit
3	Operating lamp
4	Dentist element
5	Foot controller „I-Panel“
6	Patient chair

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General

Foreword 3
 overview 4
 Table of content 6
 Technical data 8
 replaceable parts 9
 Signs and symbols 10
 Intended use/General guidelines 13
 Working area 14

Dentist element

Description 18
 Unit startup 19
 Control elements 20
 User level 24
 Air turbine 25
 Micromotor 26
 Endo-function for micromotors 29
 Ultrasonic scaler 32
 3/6-Function-syringe 36
 Camera 37
 Regulation of spray-water 38
 Second instrument status (2IS) 39
 Operating lamp - MoonLite 40
 Operating lamp - SunLite 42
 Operating lamp program 44
 Instrument rinse program 45
 Timer function 46

Patient chair

Description - Patient chair 48
 Working area 50
 Initial operation 53
 Control elements, manual control 54
 Program control 55
 Safety switches 57
 Headrest 58

Assistant element

Description assistant element 62
 Cup filler 64
 Water regulation 66
 Big suction, small suction 67
 3-function syringe (optional) 68
 Curing light (optional) 69

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Amalgamabscheider MST 1 (Optional)	72
Mundspülbeckenventil	73
Aqua-Bottle-System.....	74

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Cleaning and Maintenance

Daily	84
Quarterly	99
When required	100
Cleaning / maintenance - MoonLite	102
Cleaning / maintenance - unit	104
Readjustment of the backrest (Guidance skid)	105
Change of suction hoses	108
Fuses.....	109
Debugging	111

Mounting

General mounting.....	116
Unit transport.....	117
Patient chair installation	118
Unit column opening	122
Unit configuration	123
Overview of Appendixes.....	125

Appendix

Setting up drawing	128
Valve block.....	136
Air-water-schematic	137
Explanation of signs	140
Mounting instructions accessories	141
Safety inspection	142
EMC-declaration	144
Product information – Ritron	148
Product information – Ritren	150
Guarantee	152

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Dental work place

Typ:	Contact Lite Comfort Contact World
Rated voltage:	230 V~ 50/60 Hz
Rated input:	2300 VA
Mode of operation:	DAB*
Protection class:	I
Type of applied part:	BF
Fuse:	10 A time lag
Electrical protection: (medical rooms)	Power supply according to DIN VDE 0100-710
Water pressure (input):	3.5...5.0 bar**
Requirements for water quality: optimum of water hardness:	1.42 ... 1.78 mmol/l (8 ... 10°dh)
Air pressure:	5.5...7.0 bar
(in case of higher pressures, install a pressure reducer) (Water pressue/Air pressure configuration with pressure regulator → see Air- and water plan (U10056, U10057) in appendix of this manual)	
Max. water consumption:	
Dentist unit:	0.15 l/min
Assistant unit:	4 l/min

Amalgam separator (optional)

Typ:	Metasys Multi System Type 1
Separation efficiency factor:	at least 95 %
Discharge rate max.:	4.5 l/min
Check character:	PA-II 3816
(Institute for building technology, Berlin, 27. 11. 89)	

Patient chair

Rated voltage:	230 V~ 50/60 Hz
Rated input:	200 VA
Mode of operation:	DAB* (25 s/400 s)
Protection class:	I
Type of applied part:	B
Maximum stroke:	400 mm Lite 500 mm World
max. Load capacity	160 kg

Dental operating light „MoonLite“

Typ:	LED
Rated voltage:	17-24 V ~ 50/60 Hz
Rated input:	26 VA
Mode of operation:	continuous operation
Protection class:	II
Illuminance	3000-50000 lx
Colour temperature:	5000 K
Colour rendering index (CRI):	>97

Dental operating light „SunLite“ (optional)

Typ:	LED
Rated voltage:	17-24 V ~ 50/60 Hz
Rated input:	20 VA
Mode of operation:	continuous operation
Protection class:	II
Illuminance	8000-35000 lx
Colour temperature:	4200-5800 K
Colour rendering index (CRI):	>85
Lifetime cycle:	50000 h

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Syringe

“Ritter TopJet”

Type of applied part: BF
Mode of operation: DAB (10 s/20 s)

Syringe

“Luzzani 6F Minilight”

Type of applied part: BF
Mode of operation: DAB (10 s/20 s)

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List of removable/interchangeable parts

<u>Part</u>	<u>Description see page</u>
Filter of suction system	5
Hand piece of the syringe	5
Ultrasonic scaler	5
Suction hose	5

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Operating conditions

Temperature: +15°C to +40°C
 rel. humidity: 20% to 70%
 Barometric pressure: 860 hPa to 1060 hPa

Transport and storage conditions



(original packing)

Temperature: -5°C to +60°C
 rel. humidity: ≤90%
 Barometric pressure: 700 hPa to 1060 hPa

Icon	Explanations to type plates
	Serial number
IPX1	Closed equipment Protected against spilling water (foot controller)
IP20	dental unit (Protection against foreign bodies, no protection against water)
	CE-mark according to Directive 93/42/EEC for medical products
	The symbol of the crossed waste container indi- cates, that Ritter products are electrical equip- ment according to Directive 2002/96/EC. Within the EU-Countries these products have to be separately disposed. Ritter wants to encourage you, to minimize possible effects of waste to environment and health outside of the EU, too. Please follow the local directives for disposal of waste and make use of the separate collection of electrical equip- ment.
	Pay attention the safety instructions for this product.
	Type BF applied part
	Following the instruction for use
	Manufacturer/Date of manufacturing
	Medical device

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Our equipment complies with the requirements of standards for medical products (MPDG) and Directive 93/42/EEC for Medical Devices, revised by EU-Directive 2007/47/EG from 2010-03-21.
 It also complies with the VDE-requirements for electromedical devices DIN EN 60601, part 1 / DIN VDE 0750, part 1. It is suited for connection to supply lines installed according to the DIN VDE 0100-710 regulations for medically used rooms.

Intended use

	This dental treatment unit is used for diagnosis, therapy and dental treatment in humans by qualified personnel.
	The device is not suitable for operation in potentially explosive areas.
	This dental treatment unit must not be modified without permission of the manufacturer.
	We can only accept responsibility for the safety properties of our equipment, if installation, maintenance, repair and modifications are carried out by our workshop or by a person explicitly authorized by us. Further, the equipment has to be correctly handled according to our operating manual.
	In case of failure, structural parts, which can effect the security of the device are to be replaced by original parts only.
	In case of repair, the technician must issue a certificate as to the kind and extent of work performed and, if applicable, indicate any changes of nominal data and operating range. The certificate must contain the date of repair, the name of the executing company and a signature (according to „MPDG“, DIN VDE 0750 / DIN VDE 0751).
	Close to the equipment do not operate any devices producing electro-magnetic emissions (such as wireless phones, micro-wave therapeutics devices,...)
	When moving the back rest, the arms of the patient must not be placed behind his back! Otherwise there is a risk of squashing!

Disposal

According to law, electric equipment is to be disposed by authorized disposal firms only.

Technical controls within the scope of maintenance

Scope and schedules of safety-relevant inspections are as follows:

Every two years the following safety-relevant inspections according to DIN EN 62353 have to be performed:
 see Form „Safety inspections – Dental treatment unit“.

We recommend to keep an equipment control log for documenting the results of the safety-relevant inspections.

The safety standards for medical products and all legal obligations are to be fully observed by the user.

The unit lifetime cycle is: 10 years.

	Operator (intended user) Trained and instructed dental professionals
---	---

Biocompatibility is demonstrated according to EN ISO 10993-5:2009, EN ISO 10993-10:2010

When imported into countries whose native language is not English or German, the operating manual needs to be translated into local languages. In this translation, replace the words "original operating manual" by the following notice: "translation of the original operating manual."

The Ritter® Concept GmbH is not directly or indirectly liable for mistakes or failures in this documentation as well for the given benefits of the use with this operating instruction!

The illustrations in this manual may differ from the real unit.

The duplicating of this documentation without written permission of Ritter® Concept GmbH is prohibited!

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! Multimedia
 For components like Intraoral camera, monitor, observe the separate User manuals for these devices!

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! Only accessories recommended by Ritter Concept (see following table) may be used. The use of other accessories will void both the warranty and the conformity for this dental treatment unit

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no.	part number	description
1	F000033064	Flush cover


Abb.: Arbeitsbereich "Contact Lite Comfort C"
 Fig.: Working area of "Contact Lite Comfort C"

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The assistant unit is integrated in the unit column. No cables or supply hoses are located on the floor.

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The double-arm of the assistant unit has an extremely wide working radius. The instrument holder can be brought very close to the patient or swivelled far off in order to give the assistant space.

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
Model „S“ (swinging hoses)

- Instrument carrier with swinging lever
- smooth running swivelarm (precision needle bearing), pneumatic limited
- order of instruments free selectable
- single or double tray

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Model „H“ (hanging hoses)

- instrument carrier with hanging hoses
- smooth running swivel arm (precision needle bearing), pneumatic limited
- integrated double tray
- optional acceptance of multimedia components possible

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Model „C“ (cart-system)

- mobile cart connected with supply line to the unit
- non-slip silicone plating for acceptance of two norm-trays
- big, steerable castor for good riding quality
- optional acceptance of multimedia components possible

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


Startup

The unit is activated via the main switch (1) at the patient chair pedestal. After switching-on the control lamp of the main switch lights and the LCD-Display shows the information about initialization for about 2 seconds.


2 s

After initialization the LCD-display shows for about 5 seconds the information about the actual software versions integrated in the dental work place.


5 s

Now the LCD-display shows the current user level for about 2 seconds.

2 s

After completing the start-up procedure the dental unit is now ready for operation and at the Ritter symbol appears on the LCD-display.

Daily before the first treatment hand pieces, cup filler and spittoon bowl should be each actuate for several times, in order to rinse out stale water.

See also: Rinse program at page 45.

For safety reasons, shut-off the main switch when leaving the practice!

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Graphical LCD-display (1) used for showing different information to the user:

- the adjusted rotation speed of the micromotor, performance of the turbine or the scaler (These information are only shown when the instrument is taken out of its holder)
- the adjusted times for spittoon bowl rinsing and cup filler
- displaying of defective or active safety switches during the patient chair movement

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
Abb.: Armbrakes "H" Modell
Fig.: Arm brakes "H" model

Arm brakes (H & S models):

The handle of the tableau (2) is provided with photoelectric barriers. It releases the pneumatic arm brakes at interruption.

The dentist keyboard at model implementations of the form "Hanging hoses" differs only in the arrangement of the control keys.

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 All important functions during the treatment can be accessed and controlled via the I-Panel.

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Control elements patient chair:	
If all instruments (except the syringe) are in their holders.	
1	Backrest DOWN
2	Backrest UP
3	Patient chair UP
4	Patient chair DOWN
5	Foot lever
1+2	Automatic movement to „spitting position“ actuating again moves chair into last treatment position
3+4	Automatic movement to „exit-position“

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Control elements for instruments:	
1	Reversal of rotation direction of the micromotor or switching between PARO-, ENDO- and SCALER-function - Function is only selectable at taken out, but not active instruments.
2	Chip blower: ON/OFF
1+2	Instrument light: ON/OFF
3	At taken out instrument: activates the instrument WITH spray water- or chip blower-function (only if corresponding LED at the dentist keyboard is illuminated)
4	At taken out instrument: activates the instrument WITHOUT spray water- or chip blower-function (regardless of whether the corresponding LED on the dentists element lights up)

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


The current user level is shown permanently at the upper right corner at the LCD-display. The user level 1 is not displayed.

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The current selected user level is saved when switching of the dental unit.

Function

The software of the dental unit can manage up to 4 different user levels. The following settings are saved permanently for each user level:

- all settings for the treatment instruments (turbine, micromotor and scaler)
- 4 different treatment positions for the patient chair (this provides 16 free programmable positions)
- the pre-setted time for the timer-function
- the settings for the second instrument status (2SI) for each single treatment instrument (→ see page 39)

Operation

The user level can be changed by activating the key **“special function”**.

 All treatment instruments must be in their instrument holders !

Every actuating of the key **“special function”** increases the current user level.

After user level 4 the software returns to user level 1.

The actual activated user level is shown at the LCD-display for about 2 seconds during changing procedure.

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Startup:

Take the handpiece from its holder. The LCD-display provides the following information:

- Number of the turbine (1 or 2)
- the adjusted power in percent and by a bargraph with ten steps
- the current user level (except user level 1)

The turbine is activated via the pedals of the I-panel.

Depending on the pressure of the pedal, the power of the turbine changes from a senseful minimum to the adjusted maximum power.

The maximum power value can be adjusted via the **“Plus”/“Minus”** keys at the dentist keyboard. A short actuation increases or decreases the value in one step.

If the key is actuated longer, the value of power can be adjusted in wider steps. This is only possible if the turbine is not activated.

Functions on the dentist keyboard:

Instrument light:

- switches the light at the handpiece ON or OFF

Spray water*:

- switches the spray water at the turbine ON or OFF
- for activation press the left foot pedal of the I-Panel

Chip blower*:

- switches an constant air flow at the turbine ON or OFF
- for activation press the left foot pedal of the I-Panel

^{*)} **The functions spray water and chip blower are mutually exclusive!**

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Air turbine

Second instrument status:

- switching between two different instrument settings (→ see page 39).

 The instruments name is flashing at the LCD-display, after activating the second instrument status.

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Technical data:

Micromotor „SLM-E Type 810 LED Ritter“

(brushless motor, with LED instrument light)

Motor 1:

Rotation speed-range (normal)

100 min ⁻¹ ... 1000 min ⁻¹	in steps of 100 min ⁻¹
1000 min ⁻¹ ... 10000 min ⁻¹	in steps of 500 min ⁻¹
10000 min ⁻¹ ... 40000 min ⁻¹	in steps of 1000 min ⁻¹

Rotation speed-range (Endo modus)

100 min ⁻¹ ... 1000 min ⁻¹	in steps of 50 min ⁻¹
1000 min ⁻¹ ... 2000 min ⁻¹	in steps of 100 min ⁻¹

Motor torque-range (Endo modus)

0.3 Ncm ... 1.0 Ncm	in steps of 0.1 Ncm
1.0 Ncm ... 3.0 Ncm	in steps of 0.2 Ncm

Motor 2:

Rotation speed-range (normal)

2000 min ⁻¹ ... 10000 min ⁻¹	in steps of 500 min ⁻¹
10000 min ⁻¹ ... 40000 min ⁻¹	in steps of 1000 min ⁻¹

 **For motor 2 no Endo function is available.**

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Micromotor „Infinity LED“ (brushless motor, with LED instrument light)

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 **IMPORTANT: Please read the instrument manufacturer's operating manual carefully!**

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
Startup:

Take the hand piece from its holder. The micromotor is activated via the pedals at the I-Panel.

Rotation speed-range:


The maximum rotation speed of the motor can be adjusted in 20 steps from 1 000 min⁻¹ up to 40 000 min⁻¹. The default maximum value is shown at the LCD display.

With taken instrument you can change the maximum value via the "Plus"/"Minus"-keys. The last displayed value is stored automatically.

The changing of the maximum rotation speed can be also done with active motor.

If the respective key is actuated longer, the rotation speed can be adjusted in wider steps.

The rotation speed of the motor can be changed continuously with the pedals of the I-Panel during the treatment.

Functions on the dentist keyboard:

Instrument light:

- switches the light of the handpiece ON or OFF

Spray water:

- switches the spray water at the motor ON or OFF

Chip blower:

- switches an additional lasting air flow at the micromotor ON or OFF

Reversal of rotation direction:

- switches micromotor into anti-clockwise rotation

Second instrument status:

- switches between two different instrument settings (→ see page 39)

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Functions at the I-Panel:

right pedal:

activation of the motor WITHOUT spray water or chip blower (independent from settings at the keyboard)

left pedal:

activation of the motor WITH spray water or chip blower (dependent from the settings at the keyboard)

right pushbutton:

chip blower: ON/OFF

left pushbutton:

anti-clockwise rotation: ON/OFF

both buttons simultaneously:

instrument light: ON/OFF

Brushless motors

Function:

The brushless motors can be ordered optionally when purchasing the dental unit. They feature a higher long-life cycle and are maintenance free.

Operation:

The operation of the brushless motors is the same as of the DC-motors.
(→ see page 27)

Additional it is possible to activate the **Endo-mode** by using the **“Endo”** key at the dentist keyboard.
(→ see page 29)

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
This mode is available only for the use with following motors:
→ SLM-E Type 810 LED Ritter

Description:

The Endo-mode for motors can be used for different endodontic treatment purposes. It is possible to save the maximum rotation speed and the maximum motor torque for each single motor in 8 different programs.

The ATR (automatic reverse) function effects an automatic reversal of rotation direction if reaching the set maximum motor torque.

The ATRF (automatic reverse - forward) function effects an automatic reversing of the motor rotation direction for approximately 1 second.

After this time period the motor rotates again in its primary rotation direction.

Operation


Take the desired motor from its holder or its tray.

Press briefly the Endo-key in order to switch on or switch off the Endo-function of the motor.

At the display appear the following information:

- (1) Number of the active motor
- (2) Indication of the active Endo-mode
- (3) Current Endo-program
- (4) Maximum motor rotation speed in rpm
- (5) Maximum motor torque in Ncm

The factory made pre-set values for each program are 1000 rpm and 3.0 Ncm.

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Selecting program:

Press briefly the "Star"-key in order to switch between the single programs.

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

Setting the maximum rotation speed:



Press the "Plus/Minus"-key in order to set the value of maximum rotation speed. Selectable are values between 100 rpm ... 1000 rpm. The last displayed value is saved automatically.

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Setting the maximum motor torque

Press at first the "Endo"-key longer than 2 seconds. The value for the maximum motor torque blinks at the display.

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Press the "Plus/Minus"-key to set the maximum value of motor torque.
 Selectable are values from 0.3Ncm to 3.0 Ncm.

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

After setting the desired value, press briefly the "Endo"-key in order to end the setting mode.

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
ATR (automatic reverse) /

ATRF (automatic reverse - forward):

Press the "Rotation direction"-key longer than 2 seconds in order to switch between the functions "Normal/ATR/ATRF"
 The display shows the active function.

normal function active:

ATR function active:

ATRF function active:

normal function active:

With active ATR/ATRF function a permanent signal tone sounds during the automatic anti-clockwise rotation of the motor.
 The signal tone can be deactivated via the SoftDIP switch number3/Set3.

 **The second instrument status is not available during active Endo-mode.**

With active ATR/ATRF function the anti-clockwise rotation of the motors can not be switched manually.

 **The information about the motor torque are approximate values.**

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
Technical data:

operating mode:

DAB
 (continuous operation with intermittent load. The permissible load times are adapted to the dentist working mode) 

security-class: BF

Startup:

Take the hand-piece out of the holder. The device is activated with the pedals of the I-Panel.

Intensity:

The intensity of the scaler can be adjusted in 10 steps from minimum to maximum value. The adjusted maximum value is represented at the LCD display.

With taken scaler handpiece you can change the intensity via the **“Plus”/“Minus”**-keys. The last shown value is stored automatically. The scaler can not be regulated with the I-Panel.

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ENDO/PARO-function:

This key switches the ENDO- or the PARO-function of the scaler ON and OFF

ENDO-function: ENDO-key illuminated
 display shows "ENDO"

PARO-function: ENDO-key illuminated
 display shows "PARO"

The availability of the PARO-function depends on the manufacturer of the scaler. Please read the manufacturer's operating manual.

Spray water:

- switches the spray water at the scaler ON or OFF

second instrument status:

- switching between two different instrument settings (→ see page 39)

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Functions at the I-Panel:

right pedal:

activation of the scaler WITHOUT water (independent from settings at the dentist keyboard)

left pedal:

activation of the scaler WITH water (if activated at the dentist keyboard)

left button:

switching between PARO-, ENDO- and SCALER-function

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









	<p>Pay attention the recommended power settings for the scaler with using of the different scaler tips. → Satelec: see table on page 34. → EMS: see table on page 35.</p> <p>When non-compliance of these power settings, a increased risk of breaking of the scaler tips and a encreased risk of injury to the patient exist.</p>
--	--

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



	<p>IMPORTANT: Generally the use of scaler is contraindicated for use with: Patients with a pacemaker, defibrillator or another implantable electronic device.</p> <p>Please read the instrument manufacturer's operating manual carefully!</p>
--	--

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
























Spitze Tip		Leistungs-Einstellung power setting	
ET25			30%
ET25L			30%
ET25S			30%
ET40			30%
ET40D			30%
ETBD			30%
ETPR			100%
EX1			60%
EX2			60%
EX3			60%
EXL			60%
EXR			60%
GI-1			90%
H1			10%
H2L			10%
H2R			10%
H3			10%
H4L			10%
H4R			10%

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Maximum power settings for using scaler tips (EMS):

Spitze Tip		Leistungs-Einstellung Power setting	
A			100%
B			100%
C			100%
P			100%
PS			70%
PI			70%
RT1			70%
RT2			30%
RT3			30%
H			70%
D			70%
retro Berutti			70%
RE2			70%
SD			100%
SM			100%
SB			100%
SBd			100%
SBm			100%
PF			70%
VE			100%

Spitze Tip		Leistungs-Einstellung Power setting	
E			70%
F			70%
G			70%
SP			30%
PL1			70%
PL2			70%
PL3			70%
PL4			70%
PL5			70%
DPL3			70%
HPL3			70%
K-type files			30%
ESI file			30%

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

IMPORTANT: Please read the instrument manufacturer's operating manual carefully!

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Operation of 6F syringe "Ritter TopJet":

- 1** Key for air delivery
- 2** Key for water delivery
- 1+2** Delivery of spray

The heating function can be activated via the "F1" key. The display shows the text "heating activated" and a thermometer symbol is displayed.

For technicians only:

This function is only available if the corresponding soft DIP switch has been set to on (Set: 1, switch: 8).

Functions of 6F syringe "Luzzani 6F Minilight":

- 1** Key for air delivery
- 2** Key for water delivery
- 3** Delivery of heat adjustable at the turning handle
- 1+2** Delivery of spray

The heating of the syringe is switched at the twist grip. The activate heating is displayed by an illuminated LED at the twist grip.

Operation of the 3F syringe "RitterTopflex" & "Luzzani 3F Minimate":

The operation is the same as at the 6F syringes but without the possibility of heating the media.

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	IMPORTANT: Please read the instrument manufacturer's operating manual carefully!
---	---

Functions:

The dental unit can be additionally equipped with an intra-oral camera.

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Operation:

Press the button "F2" to get into the camera mode. The display shows "camera mode active":

It is possible to generate a freeze frame by actuating one of the foot pedals of the I-panel. Actuating the pedal again quits the freeze frame. Further information can be found in the operation manual of the camera manufacturer.


Fig.: Display camera mode

 **For technicians only:**
This function is only available if the corresponding soft DIP switch has been set to on (Set: 2, switch: 8).

	Important: Use camera only with sleeves. Cleaning according the manufacturer's instruction!
---	--

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For each instrument (except the syringe) there is the possibility to regulate the spray water quantity by hand.

The regulators can be found at the backside of the keyboard (1). The number of available regulators depends on the number of instruments.

The order of regulators depends on the order of the instruments on the tableau.

Please see also pictures above!

The lowest regulator is for the instrument on the most right side of the unit, the upper regulator is for the instrument on the most left side (S-Models)

At "H" and "C" models the order of the regulators complies with the order of the instruments at the front of the tableau.

In order to regulate the water, take the according instrument from its holder and actuate it with spray water function. Now you can adjust the water quantity by turning the according regulator as illustrated in the picture below.

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Funktion:

The software comes with a function called "second instrument status". It allows the saving of an additional set of parameters for the instrument. This contains the following parameters:

- power of turbine or scaler or the rotation speed of the micromotor
- activation of spray water or the chip blower (depend on the instrument)
- reversal of rotation direction for the micromotors
- activation of the instrument light (depend on the instrument)


By pressing only one key of the dentist keyboard, it is possible to change between the different instrument settings.

Operation:

The second instrument status is activated by actuating the key „**special function**“. The instrument needs to be taken out of its holder.

The instrument's name is flashing at the LCD-display, after activating the second instrument status. All changes to the settings will now be saved automatically.

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When actuating the key „**special function**“ again, the software returns to the first instrument status. The instrument name stops flashing at the LCD-display. All changes of the instrument parameters are also saved automatically.

This provides the possibility for a fast changing of the instrument settings during a treatment by actuating the key „**special function**“.

It is possible to program the second instrument status for all four user levels. This provides 8 different saveable instrument settings for every treatment instrument.

The complete 2IS-function can be deactivated on request by a service technician.

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The dentist working place can be equipped optionally with the operating light "MoonLite".

Operating light is fixed on flexarm mounted to the unit. It is equipped with two LEDs. (Light is passed through two secondary lenses and reflected by parabolic mirrors.)

Essential features:

- comfortable, shadow less light
- easy flexible, keeps the adjusted position safely
- removable handle, sterilizable
- mechanical dimmer
- noiseless fan

optional: Proximity sensor „infrared sensor“

- Infra-red sensor enables light to be switched on or off without direct contact, thereby eliminating the possibility of cross-infection.

Control elements:

- 1 Lever
- 2 Lens
- 3 Parabolic reflector
- 4 Handle

Startup:

The operating light can be switched on alternatively with the lever (1) or the key at the dentist keyboard.

**Operating lamp with proximity sensor
Switch ON/OFF**

To switch on and off, move the hand near to the sensor up to a maximum distance of 3 cm.

IMPORTANT:

The operating light shines only if both switches are activated. For turning off, you just have to activate one of them.

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Light intensity

The required intensity can be adjusted individually by turning the dimmer/lever (1).

Reduce light intensity

Press lever (1) and keep it to the *left* side (rear view of the operating lamp) until it reaches the desired intensity.

Proximity sensor

Place the hand near the sensor until desired intensity is reached.

 If the minimum light intensity is obtained, you will hear an acoustic signal (1 beep).

Increase the light intensity

Press and hold lever (1) to the right side (rear view of the operating lamp) until it reaches the desired intensity.

Proximity sensor

Place the hand near the sensor until desired intensity is reached.

 If the maximum light intensity is obtained, you will hear an acoustic signal (2 Beeps).

To go to the minimum intensity immediately, press and release the joystick lever (1) to the front or rear. The next time the joystick (1) is pressed and released to the front or rear, it will revert back to the previously memorised intensity.

 The lever (1) must be handled delicately to prevent breakage.

 **Each time the lamp is turned on, the light intensity will be at the level memorised when it was turned off the time before.**

 **Patients at risk (e.g. with eye conditions) are not allowed to look directly into the curing beam. Use correct protection and precautions.**

 **WARNING - reversible dazzling with wrong use of LED operating lamp - A wrong use of the LED operating lamp can lead to temporary visible disturbance, therefore, when turning on or moving of lamp, never direct light field to eyes of patient, yourself and others. Always ensure a distance of approx. 70 cm between lamp and head of patients.**

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Operating light is fixed on flexarm mounted to the unit.

The clearly geometry of reflectors, in combination with each LED source is, allows to realize a luminous flux which is homogeneous, clean and without shadows.

Essential features:

- comfortable, shadow less light
- easy flexible, keeps the adjusted position safely
- removable handle, sterilizable
- dimmer
- noiseless fan
- adjustable color temperature

Control elements:

1	Handle
2	Regulator for brightness
3	Regulator for colour temperature

Startup:

The operating light can be turned on and off with the lever.

Light intensity

The required intensity can be adjusted individually by turning the dimmer (2).

Reduce light intensity

- Rotate the brightness regulator (2) *anti-clockwise* (rear view of the operating lamp) until it reaches the desired light intensity.

Increase the light intensity

- Rotate the brightness regulator (2) *clockwise* (rear view of the operating lamp) until it reaches the desired light intensity.

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Adjustable range: 4.200 K - 6.000 K

Increase the colour temperature:

→ rotate the regulator (3) anti-clockwise

Reduce the colour temperature:

→ rotate the regulator (3) clockwise

Light field

Size: 7 cm x 14 cm at a distance of 70 cm far from the source and with well defined outlines, in order to avoid risk of patient dazzling.

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	Each time the lamp is turned on, the light intensity will be at the level memorised when it was turned off the time before.
--	--

	Patients at risk (i.e. with eye conditions) are not allowed to look into the light beam.
--	---

	WARNING - reversible dazzling with wrong use of LED operating lamp - A wrong use of the LED operating lamp can lead to temporary visible disturbance, therefore, when turning on or moving of lamp, never direct light field to eyes of patient, yourself and others. Always ensure a distance of approx. 70 cm between lamp and head of patients.
--	--

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Funktion:

The operating lamp program allows an automatic activation or deactivation of the operating lamp, when arriving certain positions of the patient chair. In standard configuration the operating lamp program is always active.

For technicians only:

To switch off the operation of the operating lamp program, soft-DIP-switch 1 Set 1 must be deactivated. See manual A1180 "soft-DIP-switches"

Operation:

The operating lamp program is configured as follows:

If arriving the programmed chair positions "P1" or "P2" the lamp will always switch on automatically.

At approaching the exit- or the spitting position after key pressure the lamp is always switched off.

If the key "Spitting position" or „Exit position“ is actuated again, the chair moves into the last position and the lamp takes the previous status. This means, that the lamp switches on again, if it was "on" at the last position, or "off" if it was "off" at the last position.

For technicians only:

If the soft DIP-switch 2 Set 1 is additionally activated, the bowl rinsing will be activated automatically when arriving the spitting position.

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


Function:


The instrument rinse program cleans the instruments over an adjusted time period. Because of hygienic aspects it is recommended to do this daily before the first treatment is done.

Operation:

For the activation it's necessary that all instruments are in their holders.


Actuate the key „chip blower“ longer than 3 seconds in order to get access to the rinse program.

The display shows the preset rinse time. The flashing arrows shows, that the instruments have to be taken out of their holders. The instruments which remain at their holders, will not be rinsed.


The value of the rinse time can be adjusted in a range from 10 to 20 seconds by using the „Plus/Minus“ keys. The last displayed value will be stored.

To start the rinse program press the key „chip blower“. The decreasing rinse time is shown at the display. In addition every 10 seconds a signal sounds.

The rinse program can be prematurely stopped at every time by pressing the key „chip blower“. At the end or abort of the rinse program a long signal sounds.

The arrow at display shows, that all taken instruments have to be placed back in their holders. After completing the program can be exit by pressing the key „chip blower“.

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

>3 s


Funktion:

The timer function can be used for different purposes during a treatment.

Operation:


In order to adjust the presetted time, the key „**timer**“ has to be actuated as long as following display appears. For that purpose all instruments must be in its instrument holders.

The preferred time can now be adjusted in steps of 10 seconds from 10 to 120 seconds by pressing the „**Plus/Minus**“ keys. The displayed value is stored permanently.

The mode can be ended by a short actuation of the key „**Timer**“.

The timer is activated by a short actuation of the key „**Timer**“. If all instruments are in its holders, the decreasing time is shown at the display and the unit makes a ticking sound. 10 seconds before the time ends, a signal sounds. At the end of the time cycle a multiple signal sounds.

If any instrument is out of its holder, only the signal sounds. The timer function can be prematurely aborted by pressing the key „**timer**“.

The timer function can be deactivated on demand by a service technician.

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We recommend you to read this service manual very carefully in order to have full knowledge of all the advantages this chair is offering you.

Contact World, patient chair

The main characteristics are:

- quickset programming for four chair positions
- maintenance-free spindle motors, very reliable and extremely quiet
- great height adjustment for working in sitting or standing position
- large adjustment range of backrest, reclined and upright backrest position possible
- backrest with integrated arm support
- slim backrest for more leg space
- additional counter-movement of seat while moving the backrest (Trendelenburg-movement) increases the comfort and avoids the "undress effect"
- comfortable seat, soft and anatomically shaped
- small, round head support, enables the dentist to get even closer to the patient (standard type).
- large adjustment angle of head support for easier treatment of upper and lower jaw
- recessed backrest for head support allows better support of smaller patients
- left armrest for patient
- movable right armrest (optional)
- Operation of patient chair can be done manually by hand via dentist or assistant keyboard or by foot via I-Panel
- In order to avoid injuries, primarily during program operation, patient chair is equipped with safety switches.

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We recommend you to read this service manual very carefully in order to have full knowledge of all the advantages this chair is offering you.

Contact Lite Comfort, patient chair

The main characteristics are:

- quickset programming for four chair positions
- maintenance-free spindle motors, very reliable and extremely quiet
- great height adjustment for working in sitting or standing position
- large adjustment range of backrest, reclined and upright backrest position possible
- backrest with integrated arm support
- slim backrest for more leg space
- comfortable seat, soft and anatomically shaped
- additional counter-movement of seat while moving the backrest (Trendelenburg-movement) increases the comfort and avoids the "undress effect"
- small, round headrest, enables the dentist to get even closer to the patient (standard type).
- headrest with double pivot
- large adjustment angle of headrest for easier treatment of upper and lower jaw
- optionally: headrest with pneumatic locking
- recessed backrest for head support allows better support of smaller patients
- left armrest for patient
- movable right armrest (optional)
- operation of patient chair can be carried out manually via dentist keyboard, assistant keyboard or via I-Panel
- In order to avoid injuries, primarily during program operation, patient chair is equipped with safety switches.
- more extensive upholstery than "Contact Lite"- patient chair

Abb.: Patientenliege: Contact Lite Comfort
Fig.: Patient chair: Contact Lite Comfort

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Initial operation:

The patient chair is activated via the main switch (1) at the chair pedestal. After switching-on, the main switch pilot lamp and the LCD indicator at the dentist keyboard are illuminated.

A single long beep indicates that the unit is ready to work

ATTENTION: For safety reasons, shut-off the main switch in the chair pedestal when leaving the surgery.

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
Control elements

The patient chair can be controlled from the dentist keyboard or with the I-Panel.

Important:

If any instrument is taken from its holder (except syringe), the control with the I-Panel is not longer possible.

Manual control via dentist keyboard:

The patient chair is controlled manually by keys **“Backrest UP/DOWN”** and **“Chair UP/DOWN”** and can be brought into any desired position.

The arrows on the keys indicate the direction of the chair movement.

Manual control via I-Panel:

The pedals and buttons of the I-Panel have the following functions:

right pedal:

chair: **DOWN**

left pedal:

chair: **UP**

both pedals simultaneously:

automatic movement to exit position

right button:

backrest: **UP**

left button:

backrest: **DOWN**

both buttons simultaneously:

automatic movement to spitting position / automatic return to last position

When moving the back rest, the arms of the patient must not be placed behind his back! Otherwise there is a risk of squashing!

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Program control:

The patient chair is provided with 4 programs at each user level. With each user level 4 different chair positions can be programmed and run up automatically.

Treatment position P1

is called by key "P1" at the dentist keyboard

Treatment position P2

is called by key "P2" at the dentist keyboard

Spitting position

is called by key "spitting position" or by actuating both buttons of the I-Panel. When actuating repeated, the patient chair returns to its last treatment position.

Exit position

is called by key "exit position" or actuating both pedals of the I-Panel

Important:


Every automatic movement will be interrupted immediately, if one of the keys for chair movement is actuated at the dentist keyboard or the I-Panel.

Programming the chair positions:

The programming of the automatic chair positions can be done by using the dentist keyboard. Please proceed as follows:

- Move the patient chair manually into the desired position.
- Actuate and keep actuated one of the four program keys on the dentist keyboard longer than 2 sec till you can hear a single beep. This beep indicates, that the programming of the position was successful.

Now the automatic program position can be called by actuating the according key at the dentist or assistant tableau, as well as by the I-Panel. When actuating a chair program where no programmed position was stored before a error message will occur at the display:

Note:

The chair programs „P1“ and „P2“ can not be called via the I-Panel.

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

Abb.: Contact World

Safety switches


In order to avoid the hazard of injuries, backrest (SIS1), and the underside of the patient chair (SIS3 and SIS4) are so designed, that upon collision the chair movement stops immediately.

The safety switch in the hinged arm of the assistant unit (SIS6) will stop the chair movement, if the arm is accidentally swivelled too far into chair direction.

Function:

If a safety switch is not connected or broken, the display shows a message with the number of the corresponding safety switch and a signal sounds.

This message can be cancelled by pressing one of the keys at the dentist keyboard or at the foot controller.

If one of the safety switches is activated, the display shows a message with the corresponding safety switch.

Additionally the symbol "warning" flashes at the display and all keys at the dentist keyboard are illuminated.

This message can be cancelled by pressing one of the keys at the dentist keyboard or at the foot controller.

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Headrest

Extension and tilting of the headrest can be easily performed with one hand only. The height adjustment of the headrest is infinitely variable (self-locking).

The double joint for wide articulation can be adjusted by the knob (1). Rotating it anti-clockwise - it is unblocked, clockwise - it is blocked.

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Optionally we offer a headrest with pneumatic locking.

For moving the headrest, you have to press the button (2) shown at the upper, right picture.

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Headrest

Extension and tilting of the headrest can be easily performed with one hand only. The height adjustment of the headrest is infinitely variable (self-locking). The extension length is approx. 15 cm.

For tilting backwards, locking device (1) has to be released by pulling. Tilting forwards can be performed without pulling the locking device. The headrest can be taken off and inserted again without problems.

Optional a headrest with double pivot is possible.

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We recommend you to read this operating manual very carefully in order to have full knowledge of all the advantages this assistant unit is offering you.

Spittoon bowl is fully integrated in unit column. Instruments are on a swivel arm.

The main characteristics are:

- ceramic spittoon bowl, colour-and scratch resistant
- spittoon bowl put on, easy to remove for cleaning
- instrument holder on swivel arm holder with double articulation with an extreme wide working radius
- safety switch at swivel arm

Optionally:

- 3-part arm for instrument holder
- provided with 3/6-function syringe
- provided with curing light

Maximum swivel range of spittoon bowl: 175°


Abb.: Schwenkbereich - Speischale (optional)
 Fig.: swivel range - spittoon bowl (optional)

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
Abb.: Helferinelement
Fig.: Assistant element

- 1 Swivel arm with double articulation
- 2 Spittoon bowl
- 3 Rinsing tube
- 4 Cup filler
- 5 Mushroom
- 6 Small suction
- 7 3-function syringe for assistant
- 8 Big suction
- 9 Instrument insert for optional equipment (e.g. curing light)

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
Control elements at the assistant panel:

Key "Exit position"

Key "Spitting position/last position"

Key "Cup filler"

Key "Spittoon bowl rinsing"

Key "Operating light"

Key "Timer"

Key "Backrest UP/DOWN" and "Patient chair UP/DOWN"

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Cup filling:


The cup filling starts when pressing the key "**Cup filler**" and stops automatically after a programmed time.

Premature stop of cup filling can be released by pressing key "**Cup filler**" once more.

The recommended filling level is two thirds of the glass or plastic cup. Overflowing water will run off through the cup support plate to the drain.

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
Programming of cup filling time:

For changing the adjusted cup filling time, press and keep pressed the key "**Cup filler**" until the following message is displayed.

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Now the filling time can be adjusted by "**Plus/Minus**" keys in a range from 1..18 seconds. After adjusting the filling time, press key "**Cup filler**" to leave the programming mode.

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
Spittoon bowl rinsing:

The spittoon bowl rinsing starts when pressing key "**spittoon bowl rinsing**" and stops automatically after a programmed time. Premature stop during bowl rinsing is possible by pressing the key once more.

Programming of spittoon bowl rinsing time:

For changing the adjusted rinse time, press and keep pressed the key "**spittoon bowl rinsing**" until the following message is displayed.

Now you can increase or decrease the rinse time in a range from 10...180 seconds.

After adjusting the filling time, press key "**spittoon bowl rinsing**" to leave the programming mode.

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Regulation of water intensity:

The water-quantity of the cup filler can be regulated by adjustment knob (2), the amount of the spittoon bowl rinsing by using adjustment knob (3).

Increase the water quantity:

- Turning the adjustment knobs anti-clockwise.

Reduce the water quantity:

- Turning the adjustment knobs clockwise.

These adjustment knobs are located behind the cover door (1) in the unit column, which has a magnetic lock.

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Big suction:

When taking the big suction from its holder, the central suction system is activated. The hand piece which remains in the instrument holder should be always closed.

Design and size of the handpiece allow easy handling with one hand. It can be slid with thumb and forefinger into the desired position without torsion of the hose.


Abb. 2: Saughandstück - Spraysauger (Dürr)
 Fig. 2: Suction hand piece - big suction (Dürr)

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Small suction:

To open the handpiece simply slide the opener with the thumb in the illustrated „open“-position.

You can also regulate the suction intensity by opening the handpiece more or less.

When returning the handpiece into the holder, the suction system is deactivated.

It is recommended to close the handpiece completely anytime when returning it to its holder. This guarantees always a full suction power if using the other handpiece.

The handpiece can be equipped with commercially available suction cannulas.

For surgical cannulas a special adaptor is necessary.

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 **Suction hose and suction handpieces are expendable materials.**

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	IMPORTANT: Please read the manufacturer's operating manual carefully!
---	--

Functions:

Water, air or spray are switched by hand. Supply of the media lasts as long as the corresponding key of handle is pressed. Either air or water can be obtained by actuating the corresponding key. Actuating both keys at the same time will supply spray. Additional some syringe models offer the function of heating the media.

Operation at example "Ritter TopFlex":

- 1 Key for air delivery
- 2 Key for water delivery
- 1+2 Delivery of spray

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- Aushärtelicht "OPM"
- Aushärtelicht "Starligt" (Mectron)

IMPORTANT: Please read the manufacturer's operating manual carefully!

Operation:

Take the handpiece from its holder. The hand piece starts working by pressing the button at the instrument. After a certain time period the instrument switches off automatically. A signal sounds after 10 seconds and at the end of the time period.

The instrument will stop automatically at inadmissible heat. However, in practise you can execute enough cycles of curing in succession.

Operating hints:

Please adhere to the composite material manufacturer's operating manual as different materials and shades have specific curing characteristics. A precondition for proper curing of the composite material is that light guide and blue filter are clean (please refer to "Maintenance"). Clean very frequently with isopropyl alcohol. Avoid any contact between the polished end of the light guide and non-cured composite material. If necessary, clean with alcohol.

Do not direct the curing beam straight into the eyes!

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The (optionally) integrated amalgam separator MST1 helps you to meet all environmental protection requirements. Contamination of waste water by amalgam is drastically reduced.

The main characteristics of the amalgam separator are:

- two-step separation in short and system- adapted intervals (sedimentation and centrifugation)
- integrated dirt-resistant air/water separation
- self-cleaning centrifuge double chambers
- degree of amalgam separation: 98.6 %
- separate processing of the two different waste streams (suction, waste water)

If the amalgam collector is filled by 95% a signal beep will sound. The collector must be exchanged immediately. However, operation is still possible. To switch off the signal beep, press the „Reset“ button (1).

When the amalgam collector is 100% full, the suction flow valve does not open any more and thus operation is made impossible.

In this case or if a disruption is existent, the signal beep can't be cancelled with the „Reset“ button.

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	<p>IMPORTANT: Please pay attention to the operating manual of the manufacturer of the amalgam separator!</p> <p>In case of a disruption please call a service technician!</p>
---	---

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Optionally, the assistant unit can be provided with a spittoon bowl valve. This valve eliminates suction noises at the spittoon bowl.

In order to warrant proper function of this valve, the automatic cleaning process has to be activated at least twice a day.

The spittoon bowl valve **(1)** and the control element **(2)**, are integrated in the column of the assistant unit.
For activating the automatic cleaning process press key on the control element **(2)**.

	IMPORTANT: The manufacturer's operating manual has to be adhered too.
---	--


Abb.: Mundspülbeckenventil Contact World/Lite
Fig.: Spittoon bowl valve Contact World/Lite

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General directives

Manufactured under observance of: EN 13485, standards for medical products (MPDG) and the directive 93/42 EEC.

This supply unit corresponds to the norm DIN 1988-100 / DIN EN 1717 for protection of drinking water to prevent a back flow of the process water (of peripheral devices and dental workplaces) to the public water net.

	<p>IMPORTANT: The period of time in which the bottle should be used is limited to 6 months. It is necessary to replace the bottle by a new one after end of the period. Otherwise the bottle can be destroyed and damages to the dental unit and/or health might arise.</p>
---	--

 **Relevant spare part numbers:**
→ see table on page 80.

Clues for disinfection of the dental unit

With help of intended use of the Ritter® AquaBottleSystem the aquiferous systems of the Dental unit can be disinfected periodically. The water disinfection must be performed with the detergent/disinfectant names RiTron. Otherwise a proper disinfection can not be guaranteed.

The following points should be noted:

- Because of the device compatibility use only the specified disinfectant RiTron.
- During disinfection with RiTron the heater unit of the syringe "Ritter TopJet" must be deactivate (on dentist and assistant keyboard).

	<p>IMPORTANT: After every successful disinfection performing RiTron must be rinsed out of the system <u>completely</u>. Scour the water system so that the RiTron will be rinsed out the water unit. That is indicated by a loss of the blue coloration of the rinsing water. When clear water comes out from the instrument hoses and from the Cup filler, the wash cycle may be terminated - not earlier!</p>
---	--

- For an optimal disinfection the disinfectant must stay 24 h or longer in the water system (e. g. during weekend).
- It's only allowed to pour the blue RiTron in a therefore specified vessel.

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
Abb./Fig. 2: Flasche/Bottle

- It must be ensured, that there is enough RiTron in the bottle before the cleaning/disinfection start. In the case of premature emptying air would be sucked in and the system would not be filled with disinfectant completely. The system would not be filled with disinfectant completely. Only a water unit, which is filled with RiTron completely can be disinfection dependable! Furthermore RiTron can be splashed out uncontrolled squirting because of air in the water hoses (for example when activating the Rinse glass filler)
- For disinfection must be performed according to hygiene regulations in dental practise (using protective gloves, protective glasses and so on).

Filling of the bottle

The following steps describe the filling of the bottle (see figure 2):

- Set the switch at the bottle support to **Off**, screw the empty bottle out of the support and fill it with RiTron from the 5l cubitainer.
- Screw the filled bottle to the support (1). Do not screw too firmly!
- Set the switch to **On**.

1	Bottle support with bottle thread and valve lifter (installation of the bottle holder is made without of the pictured cover)
2	Toggle switch position: ON
3	Toggle switch position: OFF
4	Label with information of max. period of time of the bottle

After filling the PET bottle, screw the bottle hand-tight only into the bottle holder. When the Toggle is switched on, a pressure in the bottle builds up. This is audibly due to the incoming air.

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→ The dental unit is now ready for a disinfecting process
with RiTron.


Important safety regulation

Please read the notes!

-  The bottle is pressurized! (ca. 2.0-3.8 bar,
depends on the type of device)
-  Don't manipulate the bottle with spiky
objects!
-  Use only original bottles from Ritter®
Concept GmbH.
-  Before opening the bottle it is necessary
to check that the switch (valve lifter) is
set to **Off** (→ see page 75, Fig. 2).
-  After filling the bottle with water always
screw the bottle in the correct way into
the bottle holder!
-  The bottle must be screwed only hand-
tight to the bottle support.
-  After disinfection clean the thread of the
bottle and bottle holder regularly.
-  If the bottle has a mechanical damage,
the bottle must be replaced by a new
bottle for safety.
-  For safety reasons, the use of the bottle
is limited to a period for 6 months.
After that the old pressure bottle must be
replaced by a new one.

 In non-compliance of these safety regulation the
company Firma Ritter® Concept GmbH shall not
be liable for any resulting damages!

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
The date of expiration is indicated on the security label (sticker). (see following figure)

Example (see figure):
The stamped marking on the security label indicates that the pressure

bottle may only be used until August 2018. The pressure bottle must be replaced by a new one not later than September 2018.

**After filling the bottle with water always screw the bottle in the correct way into the bottle holder.
Note: Tighten the bottle by hand only.
Clean the bottle and bottle holder regularly!**

Note your properly replacement of the pressure bottle on the supplied cleaning plan (U000016505) of your Ritter® ABS (aqua bottle system).

The pressure bottle must not be used longer than it's indicated on the security label!


**For a new pressure bottle, please contact our customer service.
Phone: +49 7351 52925-37
e-mail: service@ritterconcept.com
Relevant spare part numbers:
→ see table on page 80.**

In the case of non-compliance with the safety instructions Ritter® Concept GmbH no liability will be assumed for any ensuing damages.

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
Performing the disinfection

RiTron should be filled into the water-bearing systems on Friday or before lasting several days of work break. Ensure that enough RiTron is in the bottle before beginning any disinfecting process. For one disinfecting process 300–750 ml RiTron is needed (depending on the type of device). Every water-bearing system must be rinsed as long as the blue disinfecting agent RiTron is leaking out.

- put all instruments on maximum flow

	<p>IMPORTANT: The heating of the „TopJet“ syringe (Dentist’s and assistant’s side) must absolutely be deactivated. It can result in a burn out of the syringe heating when the bottle system is not switched-on.</p>
--	---

- Detach the turbine, the hand pieces and the angle pieces (adhere to the instrument manufacturer’s user manual).
- Detach the water connection of the Airflow-device (if existent). Make a connection to the drain of the spittoon bowl with the hose included in the scope of delivery.
- Start with the cup filling (if connected) for about three times.
- Activate the single instruments via hand or foot controller with active water.
- If a rinsing system is connected to the suction hoses or to the separation automatic, activate every suction hose for about 60 s.
- External connected devices (e. g. Airflow, if existent at the unit) are also to be rinsed!
- Every external device must basically be rinsed! (A therefore necessary plug with hose for the rinsing of the external device connection can be ordered from the manufacturer or the dental depot).
- When the disinfection procedure has been finished, switch off the dental unit and set the switch at the bottle support to the **Off** position. The bottle will be depressurized and can be filled with fresh drinking water and 1% of RiTron before starting to work.

Example for dosage:

Bottle volume: 1.5 l ⇔ add 15 ml RiTron from the measuring cup.

When working with Ritren, wear rubber/PVC gloves and tight fitting eye protection.

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IMPORTANT: All connected water hoses must be rinsed in every case! Even if they were not used in the week before. Remaining water is the ideal precondition for increasing of germs!

Placing into operation of the dental unit

The blue disinfecting agent has to be rinsed out on Mondays or after several days lasting breaks of work. Every water-bearing system must be rinsed as long as clear water leaks out.

The following steps have to be performed:

- Fill the bottle with fresh drinking water and 1% RiTron. Set the switch at the bottle support to On position.
- Switch-on the dental unit.
- Activate the spittoon bowl rinsing and rinse the spittoon bowl for about 2 min to rinse out the stagnating water in the house main water line.
- Actuate the cup filler for about 8 times.
It is useful to empty the 8 cup fillings with the suction system to check its function.

Ensure that the hand pieces and angle pieces are detached before starting to rinse. This prevents that the nozzles are being blocked with residues.

- Rinse through every single instrument hose for about 2 min.
- Externally connected devices and external water connections have to be rinsed, too!
- When clear water is leaking out of the hoses, the turbines, hand- and angle pieces can be connected. When connecting external devices, make sure that the connectors have been disinfected to prevent the entry of new germ by the quick couplings.
- Fill the bottle again with fresh drinking water. For Disinfection of the 1.5 l water in the bottle add 15 ml RiTron (\cong 1 dosing cup).
- If foam appears on the turbine spray during the treatment, the dosage of RiTron must be decreased up to 7.5 ml per 1.5 l of water (depending on the quality of water).

The dental unit is now ready for operation.

The periodic weekly disinfection should be registered into the included hygiene plan for self control.

The water outlet of the cup filler has to be cleaned in periodic intervals.

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	Only water with RiTren or ready-made RiTren may be filled in!
---	--

Error Source

No pressure is establishing within the bottle and no disinfection agent leaks out of the instrument hoses!

- Check the correct placing of the bottle at the bottle support.
- The pressure lines have a leak → Inform a technician.
- Set the switch on the bottle support to „OFF“ position, tighten the bottle a little bit more to bottle support and switch-on again.
(Attention: Don't tighten the bottle too firmly!)

RiTren was spilled out during the disinfection process!

- Clean the dirty area with available disinfection agents.
- It has to be noticed that the dirty objects/areas are impassible again the uses disinfection agent.

During the rinsing with disinfection agent, only little amounts of RiTren leaks out of the instruments!

- The air pressure of the house connection is too less so that inside the bottle not enough pressure is built up.
(Notice the limit value!)

During the rinsing with disinfection agent, only little amounts of clear water leaks out of the instruments!

- It was forgotten to detach the instruments from the hoses before starting to rinse.
Residues (biofilm) of the hoses are blocking the instruments.

Spare parts:

Component/Spare part	Order No.
Bottle support, complete	0001321032FA
Replacement bottle 1.5 l	K001123025

Consumables

Description	Order No.
RiTren: 5 l Cubitainer (One 5 l container is included in the scope of delivery!)	U000003910
RiTren: 1 l bottle with dosing cup (one bottle is included in scope of delivery!)	U000003911

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Instructions for sample taking of the local health authority

- The samples must be extracted separately at each instrument. A mixing of the water samples does not give information about which instrument is the causing the exceeding of the limit values.
- Before sample taking all hoses (including cup filler) have to be rinsed for about 2 min to fill the hoses with fresh water.
- The water leak out on positions of the endings of the hoses (with or without plugged hand and angle pieces) and the cup filler has to be disinfected by wiping. Let the disinfectant soak long enough - as prescribed.
- Extract the necessary amount of water and transport it as fast as possible to according place (approved laboratory).
- The resulting information must reported to the local health authority within the prescribed period of time.

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



***Cleaning and Maintenance/
Debugging***

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






Conformity with the hygienic standards as well as good performance of the equipment are warranted by adhering to the recommendations for disinfection, cleaning and maintenance.

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


Cleaning and care instructions

What	How	Wherewith	Hints
After each patient			
	Operating light wipe disinfection	<ul style="list-style-type: none"> DKL proPad Disinfectant wipes 	Application time/ efficacy: 5 min: bactericidal, yeasticidal, virucidal*
	Suction suck off a cup of water with big and small suction	<ul style="list-style-type: none"> water 	cold water, at least 1/2 liter
	Instruments cleaning and disinfection-	<ul style="list-style-type: none"> water Duerr FD322 	undiluted, 1 min; Don't use agents that contain chloric substances! Please follow the instructions for correct cleaning of the instruments on the following pages.
After each patient <ul style="list-style-type: none"> before workplace preparation in the treatment area after the workplace maintenance before putting on and taking off gloves before and after non-surgical treatment after using the toilet 			
	Hand disinfection <ul style="list-style-type: none"> hold the hands/forearms separately over the spray nozzle inside the appliance Each hand should be turned over so that the sprayed substance is evenly distributed. The disinfectant must be rubbed into dry skin, for the surgical disinfection up to the elbows The hands/forearms/elbows must be kept damp during the whole application process. This means that the spraying process has to be repeated, if applicable. 	Hand-disinfection automat (disinfecting solution: Aseptoman)	Application time/ efficacy: 30 s bactericide, fungicide, tuberkulozid, levurozid, limitedly virucidal (e. g. HBV, HIV, HCV, influenza virus, rota virus) Application time: 60 s Noro virus
	Hand disinfection <ul style="list-style-type: none"> Hand disinfection by rubbing: keep hands moist for 30 s with a sufficient amount of disinfectant 	<ul style="list-style-type: none"> GREEN&CLEAN HD (Metasys) 	Application time/ efficacy: 30 s bactericide, fungicide, tuberkulozid, limitedly virucidal (z. B. HBV, HIV, HCV, Herpes-simplex Virus Typ I, rotavirus)

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*) begrenzt viruzid (HBV/HIV/HCV, BVDV, Vaccina-Viren)

What	How	Wherewith	Hints
<p>Daily</p> 	<p>Painted parts clean with sponge, surface disinfection</p>	<ul style="list-style-type: none"> DKL proPad Disinfectant wipes 	<p>Application time/ efficacy: 5 min: bactericidal, yeasticidal, virucidal*</p> <p>Do not use disposable latex gloves! (Discoloration of plastic parts may occur.)</p> <p>Do not use metal cleaners!</p>
	<p>clean/disinfect the upholsteries wipe upholstery, wipe disinfection</p> <ul style="list-style-type: none"> headrest upholstery backrest upholstery seat upholstery <p>Clean/disinfect the armrest wipe the armrest, wipe disinfection</p>	<p><u>Cleaning:</u></p> <ul style="list-style-type: none"> sponge, lukewarm water, soap DKL Dental care set <p>Clean heavily soiled areas with a soft brush. After being cleaned, the upholstery should be wiped dry with a cloth.</p> <p><u>Disinfection:</u></p> <ul style="list-style-type: none"> DKL proPad Disinfectant wipes <p>We recommend applying disinfectant alcohol-free spray foam for cleaning and disinfecting imitation leather upholstery to your upholstery.</p>	<p>Remove soapy water with a damp sponge and rub dry with a soft, fiber-free cloth.</p> <p>Application time/ efficacy: 5 min: bactericidal, yeasticidal, virucidal*</p>
	<p>Dentist element, Tray, Instrument holder wipe disinfection</p> <p>Assistant keyboard, instrument holder wipe disinfection</p>	<p><u>Disinfection:</u></p> <ul style="list-style-type: none"> DKL proPad Disinfectant wipes 	<p>Application time/ efficacy: 5 min: bactericidal, yeasticidal, virucidal*</p>

 After cleaning and disinfecting the surface we recommend wiping with a clean cloth the surface to remove any residues.

*) limitedly virucidal (HBV/HIV/HCV, BVDV, vaccinia-viruses)

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Burr instruments:

- Every day before the first treatment, let instruments run for approx. 120 seconds to purge stale water at all spray tubes.
- Clean and maintain turbine, handpieces and contra-angle attachments.
- Read the chapter „Rinse cycle“ on page 45 for more information.

	IMPORTANT: Please read the manufacturer's operating manual for the burr instruments very carefully.
---	--

	IMPORTANT: Instrument hoses may not be cleaned in the thermodisinfector. Temperature load for the tubes is max. 40°C
---	---

Suction handpieces:

Detach handpiece from hose. Clean accessible hose ends and handpieces. Lubricate O-rings with vaseline and reassemble handpiece to the hose. Unsuitable and not allowed are cleaning agents based on quaternary ammonium compositions, cleaning agents containing phenol and aldehyde as well as foam generating agents.

	IMPORTANT: For dosage directions, please refer to the recommendations given for the individual cleaning and disinfection agents.
---	---


Abb.: Saughandstück - Sekretsauger (Dürr)
Fig.: Suction hand piece - small suction (Dürr)



Abb.: Saughandstück - Spraysauger (Dürr)
Fig.: Suction hand piece - big suction (Dürr)

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Flush cover - Cleaning/Disinfection of hoses

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1	for air turbine and micromotor
2	for scaler handpiece
3	for Ritter TopJet/Topflex syringe
4	for Luzzani 3F/6F syringe
5	for stopper of SLM-E Type 810 LED Ritter motor
6	for stopper of Infinity motor
7	for stopper of NSK motor
8	Flush cover

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Picture	Description	Order No.
	Flush cover	F000033064

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


Picture	Description	Order No.
	stopper for NSK motor	F000033859
	stopper for SLM-E Type 810 LED Ritter motor	F000033863
	stopper for Infinity motor	F000033861

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
Picture	Description	Order No.
	Syringe ring for Topflex/Topjet	F000033086

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Picture	Description	Order No.
	Syringe ring for Minimate	F000033087

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Picture	Description	Order No.
	Syringe ring for Minilight	F000033088

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Information about Instrument rinse program:
→ see page 45

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
Abb. 2: Position des Filters
Fig. 2: Position of the filter

Filter of suction system

ATTENTION: The amalgam rests should not be rinsed into the urban waste water!

Please use a special container for collecting the amalgam rests. We also recommend the use of protective gloves.

Procedure - Contact models:

- Open the side door of the dental unit (fig.1)
- Pull out the filter element (fig. 2,3)
- Place the filter element over a suitable container and rinse it with water
- Insert the cleaned filter element
- Close the side door of the dental unit

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
Abb. 5: Siebschublade
Fig. 5: Filter chamber

Change/clean the filter chamber

	ATTENTION: The amalgam rests should not be rinsed into the urban waste water!
---	--

Please use a special container for collecting the amalgam rests. We also recommend the use of protective gloves.

Procedure to changing/cleaning the chamber from place selection system

- open the side door (1)
- pull off the cover (2) during a slight rotatory motion
- pull off the filter chamber (3) from place selection system (4)
- place the filter chamber (3) over a suitable container and rinse it with clean water
- put the cleaned filter chamber (3) into the place selection system (4)
- slip the cover (2) on place selection system (4)
- close the side door (1)

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
Cleaning of the spittoon bowl

Before proceeding with the following steps it is absolutely necessary to switch **OFF** the unit via the main switch in the patient chair pedestal!

Remove and clean insert (1) and gold trap (2). For cleaning, the rinsing tube (3) and cup filler (4) can be removed. First remove cup filler (4), then the rinsing tube (3) can be removed together with the spittoon bowl (6). Finally, remove and clean the cup support (5).

Please do not use scouring powder, aggressive or foam generating cleaning agents.

Cleaning personnel have to be instructed not to use any common household cleaning agents for cleaning the spittoon bowl.

It is prohibited to drain chemicals (i. e., disinfectant agents, household cleaners) via the spittoon bowl.

The spittoon bowl (6) can be cleaned in a thermodesinfektor with a temperature up to 95°C.

Please lubricate O-rings with vaseline before inserting cup filler and rinsing tube.

IMPORTANT: Amalgam residues from the gold trap should not be poured into the waste water.

Cup filler:

Daily, before the first patient treatment, actuate the cup filler three times to have fresh water.

Glass cups

Glass cups have to be reprocessed according to the manufacturers guidelines.

IMPORTANT: Before taking off the cup filler and the rinsing tube switch off the unit at the main switch in the patient chair pedestal.

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
	IMPORTANT: With not of "Ritter" produced syringes please adhere to the manufacturers operating manual!
---	---

Preparation instruction for TopJet/TopFlex syringe to be performed after each patient

(according to DIN ENISO 17664)

 Warning	<ul style="list-style-type: none"> • Sterilization must be performed even when using the syringe for the first time • Inappropriate sterilization is dangerous for patients and operators • Do not perform spray sterilization • Do not perform cold chemical or hot air sterilization • The appointed staff must be specialized and trained • Use the disinfectant in accordance with the specifications set by the manufacturer indicated on the label/ instructions for use • Do not use chlorine-based solutions • When sterilizing in an autoclave, do not exceed the maximum permissible filling quantity
Preparation / pre cleaning	<ul style="list-style-type: none"> • Internal cleaning: Activate water and air of the syringe for 30 seconds to clean the internal channels • Pull the grip sleeve (2) off the syringe body. All subsequent pre-cleaning steps are to be applied to the grip sleeve (2) and the cannula (1) only. The syringe body remains on the treatment unit. • External cleaning: Remove all visible surface dirt from the handle sleeve and cannula by brushing (with a soft brush) under running water (min. drinking water quality, <38°C) for at least 10 seconds. • Remove cannula (1) from handle sleeve (2) • Carefully wipe off any contamination on the tip of the cannula with a soft, damp cotton swab or a soft, damp cloth. • Then clean the cannula and handle sleeve with a suitable cleaning agent according to the manufacturer's instructions (e.g. DÜRR ID 215, 2 % solution). After pre-cleaning, rinse thoroughly with drinking water quality water (<38°C, min. 15 seconds). • Transport the cannula and handle sleeve in a suitable container for automatic cleaning/disinfection.

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<p>automatic cleaning & disinfection</p>	<p>IMPORTANT: Use suitable adapters for cleaning and disinfection (e.g. Miele A803, A833 for the handle sleeve and ADS2 for the cannula).</p> <ul style="list-style-type: none"> Thermal disinfection according to ISO 15883-1 (90°C and 5 minutes holding time; $A_0 = 3000$; e.g. using Miele PG 8581; cleaner: Dr. Weigert neodisher MediZym) the following process steps must be observed: <ul style="list-style-type: none"> Pre-rinse: 2 minutes with cold tap water (16°C±2°C) Cleaning: Clean at 45°C with tap water and 0.5% cleaning solution for 10 minutes Intermediate rinsing: 2 minutes with cold tap water (16°C±2°C) Rinsing/disinfection: At min. 90°C for 5 minutes with demineralized water Drying: Dry at max. 90°C for 15 minutes <p>The syringe body is not suitable for this purpose!</p>
<p>automatic drying</p>	<p>Automatic drying is normally part of the cleaning/disinfection cycle of the thermo washer disinfectant and should not exceed 90 °C.</p>
<p>manual drying</p>	<ul style="list-style-type: none"> Dry the outside of the instrument with a lint-free towel. Dry the lumen of the instrument with filtered compressed air (max. 3 bar).
<p>Maintenance and checking</p>	<ul style="list-style-type: none"> no particular maintenance is required Inspect for damage or signs of wear and replace parts if necessary
<p>Packaging</p>	<ul style="list-style-type: none"> only sterilize and store in suitable packaging according to DIN EN ISO 11607-1
	<p>ATTENTION! The syringe body is not suitable for sterilization.</p> <ul style="list-style-type: none"> sterilize in class-B steam sterilizer (according EN 13060, ISO 17665-1) sterilize with 3 phases fractionated vacuum 134°C±1°C in a device according DIN EN 285 134°C ± 1°C 4 minutes holding time 20 minutes drying time the autoclave must be qualified

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Preservation	<ul style="list-style-type: none"> • no particular requirements • keeping them in their sealed and sterilized package • protect from direct sunlight • keep them in an appropriate dry environment • the storage place should be the bio-burden as low as possible
---------------------	---

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 **IMPORTANT: With not of "Ritter" produced scalers please adhere to the manufacturers operating manual !**

Ultrasonic scaler:

- Scaler tip must be replaced after each patient
- cleaning and disinfection after each patient according to the separate instruction of the scaler
- sterilize the scaler handpiece according to the separate instructions of the scaler

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Aqua-Bottle-System

Daily

1. Cleaning and disinfection of bottle holder Remove the ABS bottle in use from the bottle holder. Thoroughly clean and disinfect the bottle holder and standpipe after use with a suitable cellulose cloth soaked in undiluted BC-San 100.

2. Cleaning and disinfection of ABS-bottle For the preparation of the inner surfaces, the contents of the bottle of BC-San 100 must be filled into the emptied ABS-bottle. Afterwards, the ABS bottles are to be filled with drinking water up to the screw connection (thread / quick release)

3. Exposure time Leave the solution in the ABS bottles for at least 6 hours (overnight, for example) .During disinfection, keep the filled ABS bottles closed in the preparation room or attach them as usual to the ABS holder.

Caution: During disinfection of the ABS bottle on the unit, the ABS bottle must not be pressurized! • After the reaction time or before recommissioning, the solution with BC-San 100 should be disposed in the drain of the washbasin. • Then rinse the ABS bottle thoroughly with water. See also Instructions for Use: BC-San 100 from ALPRO Medical GmbH“

Weekly or before several days of work breaks RiTron must be initiated to flush all water-bearing paths with this disinfectant.

For more information about the „Aqua-bottle system“:
→ see page 74ff

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To be performed by service technician only!

	ATTENTION: Before any repair and replacement, shut off water line and make it non-pressurized!
---	---

Check the filter insert of air and water inlet group quarterly. In case of contamination it has to be replaced by a new filter insert.

- Filter replacement water inlet group, refer to figure (2)
- Filter replacement air inlet group, refer to figure (3)

Furthermore we recommend checking and cleaning of the filter within angle valve.

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1	regulator (air) for micromotor, turbine, bottle, and i. g. spittoon bowl valve
2	air outlet for Spittoon bowl valve only
3	air outlet for dentist element only
4	regulator (air) for syringe only

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- Empty turbine exhaust air oil container (1) and replace cotton wool. The container is located at the bottom side of the dentist unit

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- The specified obligatory protective filter (2) at site must be subjected to regular maintenance.

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	Please adhere to the manufacturer's regulations.
--	---

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Aqua Bottle System und Aqua Hygiene System:

Aqua Bottle System and Aqua Hygiene System:

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	<p>An exchange of the bottle must be performed after a mechanical destruction of the bottle or after a half year of operating time!</p> <p>Aqua Hygiene System: After one year the filter should be replaced to guarantee the operating safety of the waterways.</p>
--	---

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

	The maintenance instructions of the water sterilisation system can be found at the enclosed manual of the manufacturer.
--	--

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Swinging lever – Adjust spring force

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
- pull down the swinging lever from position (A) to the lowest position (B)
- adjust the spring force by adjusting screw (C):
 - rotate the adjusting screw clockwise
 - > swinging lever will be move more tightly
 - rotate the adjusting screw anti-clockwise
 - > swinging lever will be move smoothly


Rotate the adjusting screw by hand only with little force!

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
Swinging lever – Adjust end position

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- swinging lever in not adjusted position (D)
- adjust the end positions all of the swinging levers separat by adjusting screws (E):
 - rotate the adjusting screw clockwise
 - > end position will be move up
 - rotate the adjusting screw anti-clockwise
 - > end position will be move down
- adjust all swinging lever to the same end position (F)

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Operating light MoonLite

	It is strictly forbidden to use abrasive substances, detergents with trichloroethylene, benzene, oil of turpentine or similar on any components of the lamp.
---	---

	Failure to comply with this indication might cause: Risk of breakage of plastic parts.
---	---

Cleaning the parabolic reflectors (3)

- clean with a soft cloth and isopropyl alcohol - or using the detergent “Faro Perflex”.

	Do not use detergents containing surfactants or water repellent agents, whose build up can leave halos.
	PLEASE NOTE: other products might damage the parabolic reflectors.

Cleaning plastic parts

For cleaning plastic parts of the operating lamp MoonLite, do not use detergents or disinfectants containing:

- Ammonium Hydroxide
- Sodium Hydroxide
- Methylene Chloride
- Methyl Alcohol

following disinfectants are suitable:

- Eco Jet-1 (Cattani Group)
- Sporekin Plus DS (Ims srl)
- Zefirol Quick (Molteni Dental)

After disinfection using the product “Faro Perflex” for removing residue and margins.

Disinfection

Disinfect surfaces with alcoholic disinfectants.

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Sterilizing the handles

Handles are not provided sterile and must, therefore, be sterilised before use.

To remove the handle, unscrew button (A) and slip it off.

To insert it, push firmly and screw on (A)

Handles can be sterilized up to 200 times:
Sterilize with standard cycles at 121°C-134°C.

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Cleaning of the dental work place:

Clean painted and plastic parts with water and soap.
For disinfection and cleaning please only use agents usual in the market.

	We would like to emphasize that the manufacturer's operating manual must be strictly adhered to when handling disinfectants.
---	---

Otherwise, discoloring of covers cannot be excluded entirely.

When wearing disposable gloves made of latex, discoloring of the plastic parts may occur. No such discoloring was registered with disposable gloves made of silicone, vinyl or polyethylene.

Do not use any metal polish!
Please remove drops of water, blood or chemicals immediately!

Cleaning of upholstery and covers:

Use only warm water and soap. Remove soap residues with a damp sponge. Dry with a lint-free soft cloth.

weekly

Disinfection of the water-bearing systems

Because of the device compatibility use only the specified disinfectant RiTron.

► Information about Aqua-Bottle-System:

 **see page 74ff of this user manual**

The periodic weekly disinfection should be registered in the included hygiene plan (U16505) for self control.

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The World patient's chair was manufactured with highest quality standards and delivered to the customer.

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Nevertheless the clearance between backrest and underside of the patient's chair should be checked by a service technician from time to time (Figure 1).

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If the clearance is too high, it can result in damages to the varnish in this area. (Figure 2)

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To regulate this, the following operations have to be performed by the service technician.

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- Remove the protection cover by loosening from adhesive tape
- Loosen carefully the upholstery from the underside of the patient's chair
- Unscrew the 4 screws of the motor covering (Figure 3)
- Remove the motor covering

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- Move the backrest completely back
- Loosen the M6 counter nuts on both sides of the guidance carriage (Figure 4)

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- Readjust the backrest by tighten up slightly the socket head screws (Figure 5)

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- After effective adjustment, the backrest has only a minimal clearance

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	<p>Do not tighten up the hexagon screws too firmly, otherwise it may result in damages to the track of the guidance carriage when moving the backrest!</p> <p>If the safety switch SIS1 is permanently activated, the adjustment screw at the safety switch has to be regulated. (Figure 6)</p>
--	---

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- After completing the adjustment do not forget to tighten up the counter nuts (Figure 7)

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
Change of suction hoses (Big suction/small suction)

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
- disconnect suction hose (big/small suction)
> Pull on gray end of piece (2) only!

Installation takes place in reverse order.

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
Replacing the fuses on the mainboard:

- remove cover from patient chair pedestal
- remove cover plate

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
Replacing the fuses on the board (of patient chair):

- remove cover from patient chair pedestal

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



You can find more instructions for maintenance in our training materials.

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Message	Reason	Troubleshooting
Welcome Initialising Display...	normal startup message during the initialization phase of the dental unit	informative only, no activities required
Welcome Display Rev 1.1 Tableau Rev 1.0	normal startup message during the initialization phase of the dental unit with additional display of the installed firmware version	informative only, no activities required
user level: 2 activ!	Shows message about changing the user level (→ see page 24) possible values: 1, 2, 3, 4	informative only, no activities required
button F2 not occupied!	message about not used function key possible values: F1, F2, F3, F4	informative only, no activities required
position not programed!	the selected position for the patient chair is empty (this position is not programmed)	reprogramming the wanted position (→ see page 97)
safety- switch activ! + 2	A safety switch was activated. This can happen if objects are under or behind the patient chair or the position of the assistant tableau is too low to the patient chair. (→ see page 57) possible values: 1, 2, 3, 4, 6	<ul style="list-style-type: none"> remove any objects in the movement area of the patient chair swing out the assistant tableau press any key to clear the error message
safety- switch broken! + 3	A safety switch is defect or is not connected. possible values: 1, 2, 3, 4, 6	<ul style="list-style-type: none"> press any key to clear the error message until the error is corrected, the patient chair can be controlled manually only call service technician
footcontroller malfunction! press * button	A control element of foot switch is when turning on as defective element detected. The defective element is shown on display. (in this picture: top right button)	<ul style="list-style-type: none"> receipt the error message with "*" -key the defect control element is deactivated all other functions can be used normally call service technician

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General mounting

Installations at site such as electricity, fresh water, waste water, air etc. must correspond to the local regulations.

The connection lines at site must be installed according to the actual installation plans for the units:

Contact World U11067
Contact Lite Comfort: U11068

 **If you have these plans on a scale 1:1 they can be used as a pattern for drilling the dowel holes.**

At site there has to be an easy accessible manual shut-off valve and back-flush protective filter (80-120 µm) in the water inlet pipe to the power supply unit.

Furthermore the installation must be according to the official standard EN 1717.

Before connecting the supply unit, purge the water line and the air line at site to prevent dirt particles coming into the unit.

Electrical connection:

The supply unit for the dental work place is prepared for fixed connection. For fixed connection, the neutral wire is not to be fused.

	Attention! High voltage! Only qualified technicians are allowed to perform installation!
---	---

 **For installation of the dental treatment unit please take note of the Installation manual (Doc. No. C20-0030).**

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Transport of the water unit:

To avoid damage during the transport, it is recommended to use special carrying equipment.

The picture below shows the mounting of the carrying equipment.

At first detach the protective cover (1).

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Trage (1 Paar) (5) → Bestell-Nr.: U5991

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Transport of the patient chair

In order to prevent transportation damages, it is recommended to use a special carrying device.


How to assemble the carrying device to the patient chair is illustrated in picture below.

Before assembling the carrying device, please unscrew the patient chair from the transportation pallet (1):

Loosen and remove the two hexagon screws M10 (2), wrench size 17. Remove motor housing (3). For this purpose, remove the two pan head screws with a cross-slot (4) and pull out motor housing to the front.

Carrying device (1 pair) (5) → order no.: U5991

Adapter for carrying device (2x) (6) → order no.: U5988

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




Installation of the patient chair

1) Remove the upper and side parts of the packing while leaving the chair assembled on the pallet (1) in order that it can be transported more easily. Remove the protections of the chair and all the fittings contained in the packing itself.

2) Unscrew the nuts (2) that are fastening the ground plate to the pallet and place the chair at the point pre-arranged for the installing.

3) Place either the supplied steel disks (7) or the metal sheet in form of the ground plate of the chair (not shown at the picture!) under the levelling screws (6) of the ground plate.

After this, the ground plate of the chair (5) can be mounted to the floor. Therefore, according to the kind of material of which the floor is made (lumber, concrete), correlative dowels and screws (2x) with a minimum diameter of 10mm are needed. Before the final fastening of the ground plate (5) on the floor, the level of the chair can be adjusted by using the levelling screws (6).

If necessary, act on the rear screws (3) and remove the guard (4), unscrewing the screws and pulling it upwards.

Assembly of left armrest


4) Fit the armrest on the proper bracket after having placed 2 wide washers (8). Then tighten completely the screws (10) interposing the washers (9).

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 **Attention!**
For packing reasons, the armrest is assembled downwards during the delivery of the chair!

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5) After the installation and the test of the complete unit, the front cover of the patient chair (11) can be assembled by using the allocated screw (12).

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
In order to perform the connections (for electricity, air, water) between unit column, instrument table and patient chair, several covers have to be removed when installing the dental work place. These are cover (3), lateral cover (1) and (if necessary) radius cover (5).

When removing lateral cover (1) please observe the following:

- Take off cover (3)
- Loosen cylinder screw (2).
- Open unit column door (4). Loosen cylinder screw (7) and push in holding angle (6).
- The lateral cover (1) can now be removed.

The radius cover (5) can be removed after loosening the two cylinder screws (8).

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The SoftDIP-switches are being used for the configuration of the unit. With these switches it is possible to activate or deactivate different software functions according to the customers request.


There are 24 SoftDIP-switches available, which are arranged in 3 "Sets" to respective 8 SoftDIP-switches. Depending on the installed software, not all the switches are in use.

You can get an overview about the switches and their functions from our technical hotline. Therefore we need the information about the software version, which is installed inside the unit. You can contact us also by using the following email address:
 service@ritterconcept.com

In order to enter configuration mode for SoftDIP-switches, simply press and keep pressed together "plus" and "minus"-button during you switch ON the unit with the main power switch.

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
Anzeige:

After the power-on procedure is finished, the display shows the first set with respective 8 SoftDIP-switches and the first switch is selected for modification.

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
The **“plus”** and **“minus”**-buttons serve to select the wanted SoftDIP-switch.

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With the chair control buttons the actual selected switch can be set in ON or OFF position.

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With the backrest control buttons the actual set of SoftDIP-switches can be changed.

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By pressing the **“star”**-button short, your current SoftDIP-switch setting will be saved, which is signalized by an acoustic signal, and the unit will proceed with the normal function.

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Hint:

If you are not sure about the modification you have done, you can everytime leave the configuration mode WITHOUT any modification simply by switching OFF the main power of the unit without pressing the **“star”**-button before. In this case all modifications are lost and the initial state is restored.

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IMPORTANT!

Every modification of SoftDIP-switches is without effect, unless you have switched OFF and ON the unit with the main power switch. Do NOT switch OFF the unit after you’ve pressed the **“star”**-button, unless you have heard the acoustic signal for saving the complete data, otherwise modifications may be not stored correctly.

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
	<p>AFTER you have heard the acoustic signal for saving the SoftDIP-switch settings completely, switch OFF and ON main power in order to make your modifications take effect!</p>
--	---

	Plan/Drawing	Seite/page
<u>Aufstellungspläne/Setting up drawings</u>		
<u>Contact World</u>		
Contact World S	U000012112.....	122
Contact World H/HC	U000012113.....	123
Contact World C	U000012115.....	124
<u>Contact Lite Comfort</u>		
Contact Lite Comfort S	U000012114.....	125
Contact Lite Comfort H/HC	U000012117.....	126
Contact Lite Comfort C	U000012116.....	127
<u>Installationspläne/Installation drawings</u>		
Contact World	U000011067.....	128
Contact Lite Comfort	U000011068.....	129
<u>Ventilblock/Valve block</u>		
Contact World/Lite	U000010053.....	130
<u>Luft-Wasser-Pläne/Air-water-schematics</u>		
Standard-Version	U000010052.....	131
WEK-System	U000010055.....	132
ABS-System	U000010056.....	133
<u>Montageanweisung Zubehör/Mounting instructions accessories</u>		
Contact World S+H / Contact Lite Comfort S+H	A1169	135
<u>STK-Protokoll/Safety inspection protocol</u>		
STK-Protokoll	C31-0044.....	136
Protocol for safety inspection	C31-0044.....	137

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













Acronym	Description
A	pipe socket (backwater)
AF1	solid collector 1
AF2	solid collector 2
BP1	arm brake 1
BP2	arm brake 2
DQL	compressed air source
DQW	pressure source water
DR2	throttle for spittoon bowl
DR3	throttle for cup filler
DR4	throttle for scaler
DR5	throttle for micromotor 1
DR6	throttle for turbine 1
DR7	throttle for micromotor 2/turbine 2
FLG	filter (air)
FWF1	filter (water)
GR1	supply unit A1
GR2	cuspidor
JMK	valve - cooling air
JMS	valve - spray
JMW	valve micromotor 1 (water)
JT1S	valve turbine 1 (air)
JT1W	valve turbine 1 (water)
JT2S	valve - turbine 2/micromotor 2 (air)
JT2W	valve - turbine 2/micromotor 2 (water)
JMV	valve - scaler (water)
M1	micromotor 1
M2	micromotor 2
ML	pressure gage (air)
MW	pressure gage (water)
RL1	regulator 1 (air)
RL2	regulator 2 (air)
RW1	regulator 1 (water)
SA	Small suction
SE	Big suction
SGF	cup filler
SiV1	safety valve
SP1	Topjet
SP2	Topflex syringe
SSS	sanicut bowl flush
T1	turbine 1
T2	turbine 2
VAL	shut-off valve (main inlet air)
VAW	shut-off valve (main inlet water)
VTW	manifold (water)
X1	solenoid valve (water)
X2	solenoid valve (air)
X3	solenoid valve spittoon bowl (water)
X4	solenoid valve cup filler (water)
X5	solenoid valve (air)
Y4-1	drive air
Y4-2	cooling air
Y807	valve arm brake
ZEG	scaler

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Inspection protocol for safety inspection according DIN EN 62353			
Client (dental practice)	Contractor		
Device type:	Manufacturer:		
ID No.:	Insulation class: 1		
Serial No.:	Type of applied parts: BF		
Safety inspection according: EN 62353/VDE 0751	Date of safety inspection:		
Visual inspection			
<input type="checkbox"/> Protective earth conductor not damaged (for insulation class 1 only)			
<input type="checkbox"/> Cover and mechanical parts not damaged			
<input type="checkbox"/> Insulating parts not damaged			
<input type="checkbox"/> Device connecting cables including plugs without defects			
<input type="checkbox"/> Labels exist or completed, nameplate correctly			
<input type="checkbox"/> Others			
Measurement			
Measurement	measured value	limit value	OK
Protective earth resistance		<0.300 Ω	
Insulating resistance			
Test voltage			
Device leakage current (Equivalent device leakage current)		<10.000 mA	
Encl. leakage current (Housing leakage current)			
Leakage current of application part (Equivalent patient leakage current)			
Test voltage leakage current (nominal Voltage ±10%)		min. V max. V	
Functional testing			
<input type="checkbox"/> Function and safety testing without defects			
<input type="checkbox"/> The unit has safety deficiencies:			
<input type="checkbox"/> Short-circuit L-N			
<input type="checkbox"/> Body circuit L/N-SL			
Test cycle: months		
Next test date:		
Used measurement device	Type:		
	Model:		
Tester:			
..... Place, Date	 Signature	

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guidelines and manufacturers declaration - electromagnetic emission		
The unit is dedicated for the operation in below specified electromagnetic environment. The customer or the user of the unit must make sure, that it will be used in such an environment.		
emission-measurements	congruities	electromagnetic environment - guidelines
high frequency emission according CISPR 11	group 1	The unit uses high frequency energy exclusive for its internal function. From there the high frequency emission is very low and it is improbable, that other electrical devices are being disturbed.
high frequency emission according CISPR 11	class B	The unit is determined for the usage in all facilities including living areas and such, which are directly connected to a public power supply that supplies also buildings which are used for living purposes.
overvibration according IEC 61000-3-2	class A	
voltage fluctuation/ flicker according IEC 61000-3-3	compliant	

guidelines and manufacturers declaration - electromagnetic stability			
The unit is dedicated for the operation in below specified electromagnetic environment. The customer or the user of the unit must make sure, that it will be used in such an environment.			
stability tests	IEC 60601-test level	congruity level	electromagnetic environment - guidelines
electrostatic discharge (ESD) according IEC 61000-4-2	± 2 kV contact discharge ± 4 kV contact discharge ± 6 kV contact discharge ± 2 kV air discharge ± 4 kV air discharge ± 8 kV air discharge ± 2 kV ± 4 kV ± 6 kV indirect discharge onto vertical coupling plate (VCP)	± 2 kV contact discharge ± 4 kV contact discharge ± 6 kV contact discharge ± 2 kV air discharge ± 4 kV air discharge ± 8 kV air discharge ± 2 kV ± 4 kV ± 6 kV indirect discharge onto vertical coupling plate (VCP)	Floors should be made in wood or concrete or should be provided with ceramic tiles. If the floor is provided with a synthetic material, the relative air humidity must be at least 30%.
Radiated RF disturbance according to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.17 \cdot \sqrt{P/W}$ for 80 MHz to 800 MHz $d = 2.33 \cdot \sqrt{P/W}$ for 80 MHz to 2.5 GHz with P as the maximum of rated power of the transmitter in watt according to the transmitter manufacturer and d is the recommended separation distance in meters (m). The field strength in each frequency of stationary radio transmitters should be - in accordance with a site survey - less than the compliance level. In the vicinity of equipment that is marked with the logo „Non-ionizing radiation“, disruptions are possible.
fast transient electrical disturbance/bursts according IEC 61000-4-4+A1	± 2 kV for main cable ± 1 kV for entrance- and outgoing wires	± 2 kV for main cable ± 1 kV for entrance- and outgoing wires	The quality of the supply voltage should comply with a typical business or hospital environment.
withstanding voltage (surges) according IEC 61000-4-5	± 0.5 kV ± 1 kV Power supply port line to line ± 0.5 kV ± 1 kV ± 2 kV Power supply port line to ground	± 0.5 kV ± 1 kV Power supply port line to line ± 0.5 kV ± 1 kV ± 2 kV Power supply port line to ground	The quality of the supply voltage should comply with a typical business or hospital environment.
Conducted HF disturbance according to IEC 61000-4-6	3 V _{off} 150 kHz to 80 MHz	3 V	Portable and mobile RF communications equipment should be used in no less distance to Contact Lite Comfort including cables, than the recommended separation distance calculated from the equation for the transmission frequency equation. Recommended separation distance: $d = 1.17 \cdot \sqrt{P/W}$

guidelines and manufacturers declaration - electromagnetic stability			
voltage irruptions, short time intermittent and fluctuations of supply voltage according IEC 61000-4-11	<p>< 5% UT (> 95% irruption of UT for ½ period)</p> <p>40% UT (60% irruption of UT for 5 periods)</p> <p>70% UT (30% irruption of UT for 25 periods)</p> <p>< 5% UT (> 95% irruption of UT for 5s)</p>	<p>< 5% UT (> 95% irruption of UT for ½ period)</p> <p>40% UT (60% irruption of UT for 5 periods)</p> <p>70% UT (30% irruption of UT for 25 periods)</p> <p>< 5% UT (> 95% irruption of UT for 5s)</p>	<p>The quality of the supply voltage should comply with a typical business or hospital environment.</p> <p>If the user of the unit wants a continuously function also at the appearance of interruptions, it is advisable to supply the unit by an uninterruptable power supply (UPS) or a battery.</p>
magnetic field by supply frequency (50 Hz) according IEC 61000-4-8	3 A/m for 180s per axes	3 A/m for 180s per axes	Magnetic fields at the supply frequency should comply with the typical values in business or hospital environments.
<p>REMARK U_r means the alternating supply voltage before the exertion of the test levels</p> <p>REMARK 1: At 80 MHz and 800MHz, the higher frequency range applies.</p> <p>REMARK 2: These guidelines may not applicable in all cases. The Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> <p>Field strengths from fixed transmitters, such as base stations for radio telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast can not be predicted theoretically with accuracy.</p> <p>To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Contact Lite Comfort is used, exceeds the applicable RF compliance level above, the Contact Lite Comfort should be observed to verify intended use. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Contact Lite Comfort.</p> <p>Over the frequency range of 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Recommended separation distances between portable and mobile RF communications equipment and the Contact Lite Comfort				
The Contact Lite Comfort is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Contact Lite Comfort can thereby help prevent electromagnetic interference by observance the minimum distance between portable and mobile RF communications equipment (transmitters), depends of the output power of the communications equipment, and the Contact Lite Comfort.				
protective distance depending on the transmission frequency in m				
Rated power of transmitter P in W	150 KHz to 80 MHz	80 MHz bis 800 MHz	800 MHz bis 2.5 GHz	
	$d=1.17 \cdot \sqrt{P/W}$	$d=1.17 \cdot \sqrt{P/W}$	$d=2.33 \cdot \sqrt{P/W}$	
0,01	0.12	0.12	0.23	
0,1	0.37	0.37	0.74	
1	1.17	1.17	2.33	
10	3.70	3.70	7.37	
100	11.70	11.70	23.30	
For transmitters whose maximum rated power in the above table is not specified, the recommended separation distance d in meters (m) can be calculated using the equation of the transmitter, where P is the maximum rated power of the transmitter manufacturer.				
REMARK 1: At 80 MHz and 800MHz, the higher frequency range applies.				
REMARK 2: These guidelines may not applicable in all cases. The Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				
The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range from 80 MHz to 2.5 GHz are intended to decrease the likelihood that mobile / portable communications equipment could cause interference if it is inadvertently brought into patient areas . For this reason, the additional factor of 10/3 is used in calculating the recommended separation distances in these frequency ranges.				

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Hersteller: ALPRO MEDICAL GmbH

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RiTron

For water decontamination and for preventing the formation of Biofilm in dental units.

- **RiTron** was developed for use in “Aqua Bottle System” and in “Aqua Hygiene System” of ALPRO.
- Simple and easy to use design provides a maximum of safety and durability
- **RiTron** provides on days or overnight for sterilization and permanent removal of the biofilm in the piping (including instrument hoses and cup filler) in all dental units.
- **RiTron** pays for itself. For a disinfection cycle - on weekends or overnight - no more than 300ml-750ml are needed (depending on model). This creates a clear, predictable cost, regardless of the number of patients per day.
- **RiTron** is bactericidal, fungicidal, dissolves and removes safely and reliably. (use at least 1 time per week) both organic and inorganic residues such as lime.
After complete removal of the biofilm on the pipe walls germs have no place to hide!
- **RiTron** can long remain in the waterway, such as several weeks of Holiday time.
It thus ensures that fear reliability and sterility without the need to:
„8:00 clock in the morning - and you go home“!
- **RiTron** effective against a broad spectrum of germs and prevent rapid formation of biofilm.
The ready RiTron is extremely sparing on the material and non-corrosive!
- **RiTron** was Mikrobiol its microbicidal effectiveness of the Institute for. Investigations W. Niedermann, Überlingen (2014-05) examined.

Packing unit: CubiTainer 5l **RiTron**
Order/Part No: U000003910, capacity: 5l

We look forward to your order:
Ritter®Concept GmbH
Freiburger Str. 45 · 88400 · Biberach
Germany

Fon: +49 7351 52 925-37
Fax: +49 7351 52 925-21
info@ritterconcept.com
www.ritterconcept.com

Manufacturer: ALPRO MEDICAL GmbH

Technical modifications reserved!

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Hersteller: ALPRO MEDICAL GmbH

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RiTren

For optimal hygiene of process water in dental unit water lines.

- **RiTren** was developed for use for continuous disinfection of process water in dental unit water lines (incl. instruments and cup filler) for all dental treatment units and it prevents algae formation.
- **RiTren** effective against a broad spectrum of germs
- **RiTren** prevents the formation of lime deposits and (in heating the water and longer service life) - at a water hardness of 2°dH (0.36 mmol/l) up to a water hardness of 20°dH (3.57 mmol/l).
- **RiTren** is due to the greatest purity raw materials of the highest quality, physiologically tested and classified as harmless; It conforms to the requirements of the DAB and the German Cosmetics Regulation.
- **RiTren** is non-corrosive.
- **RiTren** has a very pleasant taste (mint) and promotes the patient's well-being.
- Precondition for the application of RiTren is the Aqua-Bottle-System with 1% RiTren process water.
- **RiTren** was Mikrobiol its microbicidal effectiveness of the Institute Fresenius, Stockach 10/95 and 08/96 (RiTren) examined.

Packing unit: 1 l **RiTren**

Order/Part No: U000003911 capacity: 1 l

We look forward to your order:
Ritter®Concept GmbH
Freiburger str. 45 · 88400 · Biberach
Germany

Fon: +49 7351 52 925-37
Fax: +49 7351 52 925-21
info@ritterconcept.com
www.ritterconcept.com

Manufacturer: ALPRO MEDICAL GmbH

Technical modifications reserved!

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Terms of Guarantee

You are now owner of a high quality Ritter product. Thank you very much for your trust in Ritter Concept GmbH and we would like to maintain this! You have the warranty right of 12 months. We point to the warranty terms in our General Terms and Conditions.

We warrant, that our products are free of faults by the time of handover to the supplier. In particular, we warrant, that the production is taking place based on the acknowledged technical rules and the relevant regulations. We also warrant, that the product characteristics are available, which are determined in the contractual documents. In line with the warranty we perform a remedy of defects, when substances, parts or performances are identified as flawed or contrary to contract. The warranty starts with the date of delivery. Deficiencies are to be communicated with Ritter Concept GmbH immediately after receipt of goods and immediately after discovery. The customer is permitted to send back products without prior approval of Ritter Concept GmbH. The customer has to select the most economic way of transport.

No warranty claims exist in the following cases:

- using the product outside of or contrary to its intended use, or operating the product outside of its specifications
- failure to observe the instructions of use, maintenance or installation
- improper handling or neglected maintenance
- use of non-approved spare parts or accessories
- product modifications or conversions performed by customer or third parties
- usual wear and tear or other changes to or deterioration of the product attributable to:
 - external influences (e. g. knocks, blows, mechanical impacts)
 - water, fire
 - improper method of storage, handling or installation
 - unusual ambient conditions, especially conditions upon receiving or operating conditions at the installation location or force majeure calcification or corrosion resulting from neglected maintenance or improper handling
 - contaminants in the air and water supply
 - unusual or impermissible chemical or electrical influences
 - soiling, discolouration or wear and tear on padding

The warranty as a rule excludes parts subject to normal wear, especially lamps, glass and rubber components and cables as well as plastic parts and surfaces with regards to their colour-fastness.

All of our obligations under the warranty will be extinguished if the delivered goods are modified, improperly handled or processed. We do not warrant co-supplied or integrated third-party products. However, if requested we will assign to the customer our warranty claims against the vendor of the third-party product.

In order to activate the warranty, it is necessary to send a copy of the record of delivery and the „Warranty card“ after installation to us.