

The GIRARDELLI X-1 table darkroom meets all today's requirements of an efficient hand-developed device.

Table Darkroom

- Each part used is extremely stable and chemically resistant
- Minimal service costs over many years
- · Ideal cup size
- Allows the development of intraoral film sizes
- Even occlusal films can be developed
- Minimal maintenance
- Only short cleaning of the cups necessary when changing the used GIRARDELLI powder development concentrate
- · Simplest operation
- The film holders can be loaded and the films can be developed via the manual interventions
- Slanted top allows for excellent visibility

Technical data

Dimensions (wxdxh): Contents of the cups:

Weight:

Standard equipment:

46x32x30 cm 4x0,225 l approx. 2,3 Kg

detachable red light lens, 3 film

holders for single films,

4x0.225 litre plastic cups with

screw lid closure





Ritter worldwide -

Technicians and spare parts service!

Ritter Concept GmbH Freiburger Str. 45 88400 Biberach Germany

Fon +49 (0) 7351 / 52 925-10 Fax +49 (0) 7351 / 52 925-11

sales@ritterconcept.com www.ritterconcept.com Your Ritter Partner: